<p><b>Symbols Used on Labeling</b> Auf der Verpackung verwendete Symbole</p> <p>Symboles utilisés sur l'étiquette</p> <p>Simboli usati sull'etichetta</p> <p> Símbolos utilizados en las etiquetas</p> <p> Símbolos utilizados na embalagem</p>
---

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<b>English</b>	<b>DFU-0147-2r0</b>
<b>A. DEVICE DESCRIPTION</b> The TightRope® devices are comprised of various configurations of one or two metal buttons, a metal or bioabsorbable anchor, and suture. Single buttons (Clavicle and Coracoid) and washers may be provided in individual packages without sutures. These devices may be preloaded onto an inserter, or may be packaged with various ancillary instruments to aid in insertion. The ACL/PCL/BTB/RT/ABS TightRope® products are comprised of suture with or without a button, wedge or inserter. TightRope® ABS Implants are used with TightRope® ABS Buttons (Standard, Round and Oblong buttons) and TightRope® Button Extender. The <b>ancillary Buttress Plates are metal plates with two or four holes, depending on design, for TightRope or screw fixation.</b>	

**B. INDICATIONS**  
The TightRope devices are intended as adjuncts in fracture repair involving metaphyseal and periarticular small bone fragments where screws are not indicated, and as adjuncts in external and intramedullary fixation systems involving plates and rods, with fracture braces and casting.  
Specifically, they are intended to provide fixation during the healing process following:

- (Mini TightRope, Mini TightRope FT and CMC Mini TightRope): Hallux Valgus reconstruction (correction) by providing for the reduction of 1<sup>st</sup> metatarsal-2<sup>nd</sup> metatarsal intermetatarsal angle.
- (Mini TightRope, Mini TightRope FT and CMC Mini TightRope): Tarso-metatarsal (TMT) injury, such as fixation of foot soft tissue separations due to a Lisfranc injury (Midfoot Reconstruction).
- (Mini TightRope, Mini TightRope FT and CMC Mini TightRope): Syndesmotic trauma, such as fixation of dorsal distal radioulnar ligament (DRUL) disruptions.
- (Syndesmosis TightRope): Syndesmotic trauma, such as fixation of syndesmosis (syndesmosis disruptions) in connection with Weber B and C ankle fractures.
- (Acromioclavicular (AC) TightRope, AC GraftRope, Twin-Tail TightRope, and Dog Bone Button used with FiberTape): Syndesmotic trauma, such as acromioclavicular separation due to coracoclavicular ligament disruption.
- (Arthrex recommends the use of #5 FiberWire® or equivalent with the Clavicle and Coracoid buttons and washers.

The Arthrex Mini TightRope, Mini TightRope FT and CMC Mini TightRope, when used for fixation of bone-to-bone or soft-tissue-to-bone, are intended as fixation posts, distribution bridges, or for distributing suture tension over areas of ligament or tendon repair. Specifically, the Arthrex Mini TightRope, the Mini TightRope FT and the CMC Mini TightRope are indicated for Carpal-Metacarpal (CMC) joint arthroplasty as an adjunct in the healing process of the reconstruction of the ligament at the base of the thumb metacarpal by providing stabilization between the base of the first and second metacarpal when the trapezium has been excised due to osteoarthritis. The Mini TightRope, Mini TightRope FT and CMC Mini TightRope are also indicated for use as adjuncts in the suspension of the thumb metacarpal during the healing process of hemataoma distraction arthroplasty by providing stabilization at the base of the first and second metacarpal when the trapezium has been excised due to osteoarthritis.

The Arthrex ACL/PCL/ BTB/RT/ABS TightRope are used for fixation of bone to bone or soft tissue to bone, and are intended as fixation posts, a distribution bridge, or for distributing suture tension over areas of ligament or tendon repair. Specifically used for ACL/PCL repair and reconstruction.

- Arthrex recommends the use of #2 or #5 FiberWire® or equivalent with the Suture Button.

- C. CONTRAINDICATIONS**
- Insufficient quantity or quality of bone.
  - Blood supply limitations and previous infections, which may retard healing.
  - Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
  - Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
  - Any active infection or blood supply limitations.
  - Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
  - The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb, or disrupt the growth plate.
  - Do not use for surgeries other than those indicated.

- D. ADVERSE EFFECTS**
- Syndesmosis TightRope, Mini TightRope, Mini TightRope FT and CMC Mini TightRope Only:** Osteomyelitis surrounding the TightRope.
  - Syndesmosis TightRope, Mini TightRope, Mini TightRope FT and CMC Mini TightRope Only:** Radiastiasis resulting from a failure of the implant insertion technique.
  - Syndesmosis TightRope, Mini TightRope, Mini TightRope FT and CMC Mini TightRope Only:** FiberWire polyethylene wear-related painful aseptic osteolysis.
  - Infections, both deep and superficial.
  - Foreign body reactions.
  - Bio only:** Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
  - Silicone sensitivity, though very rare, has been reported.

- E. WARNINGS**
- The ACL/PCL/ BTB/RT/ABS TightRope devices are not intended to be used as a ligament replacement.
  - Do not add additional suture to the ACL TightRope. The extra suture may impede passage of the device through the femur.
  - Do not re-sterilize this device. Note: Metal implants which are provided alone can be re-sterilized if their sterility is compromised.
  - All metallic implant devices used for this surgical procedure should have the same metallurgical composition, including the plates that might be used in conjunction with the Arthrex TightRope Syndesmosis Repair Kit implants. Therefore, different material implants are provided.
  - AC TightRope Repair Kit, AC GraftRope, Twin Tail TightRope, Dog Bone Button with FiberTape:** Do not use this device as the sole means of reconstructing a chronic acromioclavicular joint dislocation.
  - Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device and bone.
  - Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
  - Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
  - Detailed instructions on the use and limitations of the device should be given to the patient.

- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
- Bio only:** Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.
- Once open, discard any unused suture.

- F. MRI SAFETY INFORMATION**
- This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating, migration or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.

- G. PRECAUTIONS**
- ACL/PCL/ BTB/RT/ABS TightRopes only:** Excessive force on the shortening suture strands may break the strands and impair ability to fully seat the implant. No additional force on the shortening suture strands is required when the graft/wedge construct reaches the desired position in the femoral socket and the graft stability is verified by pulling distally on the graft.
  - Load the ACL/PCL TightRope Button over the unspliced, thinner portion of the ACL/PCL TightRope suture to assist assembly.** Once assembled slide the ACL/PCL TightRope Button down to the thicker, spliced portion of the ACL/PCL TightRope suture to help prevent disassembly.
  - Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. Or, contact your Arthrex representative for an onsite demonstration. The Arthrex website also provides detailed surgical technique information and demonstrations.

- H. PACKAGING AND LABELING**
- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
  - Contact Customer Service if the package has been opened or altered.

- I. STERILIZATION**  
This device is supplied sterile. Refer to the package label for the sterilization method.

- J. MATERIAL SPECIFICATIONS**  
Refer to the package label for the materials.
- Anchor:** The anchor is made of Titanium or Poly (L-lactide acid) PLLA **Buttons/Plates:** The buttons and plates are available in two different metals in order to match the metallurgic composition of implants that may be used in the surgical procedure. The materials available are titanium and stainless steel.
- Needles (if applicable):** Stainless steel.
- Sutures:** Ultra High Molecular Weight Polyethylene (UHMWPE), or Ultra High Molecular Weight Polyethylene and Polyester (e.g. FiberWire®, TigerWire®, FiberTape®, and TigerTape sutures). All sutures meet or exceed U.S.F. and European standards for non-absorbable surgical sutures (FiberWire®, TigerWire®, FiberTape®, TigerTape suture do not meet U.S.P. diameter requirements).

- Arthrex recommends the use of #2 or #5 FiberWire® or equivalent with the Suture Button.
- Additional materials to the sutures include silicone elastomer coating (except suture with the suffix -Tape), and may include cyanocrylate and/ or nylon. The coating acts as a lubricant for suture sliding, knot tying, and ease of passing suture through tissue. Dyes may include: D&C Blue No. 6, D&C Green No. 6, and Logwood Black.**
- #2 Suture:** Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester.
- Wedge (if present):** Polyetheretherketone (PEEK).

- K. STORAGE CONDITIONS**  
Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.
- Non-bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

- L. INFORMATION**  
Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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