HATriC™

Resorbable bone Substitute

Manufactured by :

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Direction For Use HATriC™

HATriC™ is a bone void filler or bone graft substitute. HATriC™ is designed as a block of two-phase ceramic calcium phosphate consisting of 60% hydroxyapatite and 40% β -tricalcium phosphate in a wedge configuration in various sizes.

HATriC™ presents as a multidirectional interconnected porosity structure, similar to that of human cancellous bone. The HATriC™ implant slowly resorbs during the remodelling and bone defect repair process and is progressively replaced with physiological bone. The progressive resorption of HATriC™ resorbable bone void filler is intended to prevent premature resorption.

A. INDICATIONS

HATriC™ is intended for use as a bone void filler or bone graft substitute in bony voids or gaps which are not intrinsic to the stability of the bony structure (i.e. extremities, pelvis, or spine.) These bony voids and gaps may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. HATriC™ may be used as bone graft substitute for autograft or allograft when these grafts are undesirable or unavailable.

When used for posterolateral fusions, the HATriC™ must be used with autograft bone material and the ratio of autograft bone material to HATriC™ must be in the range of 1.25cm³ to 2.4cm³ of autograft bone to 1cm³ of HATriC™.

HATriC™ is a bone void filler or bone graft substitute that resorbs and is replaced with bone during the healing process. HATriC™ is not intended to be a load-bearing device. Therefore, rigid fixation techniques may often be recommended.

In addition, when used with Arthrex Opening Osteotomy Systems' or substantially equivalent systems' plates and screws, HATriC™ is intended to be used as a bone void filler or bone graft substitute in femoral or tibial osteotomies. HATriC™ is to be used in association with adequate post-operative immobilization.

B. CONTRAINDICATIONS

- Established or imperfectly dried out osteomyelitis.
- Necrosis at the recipient site.
- Degenerative bone disease in the limb to be treated.
- Intra-articular implantations.
- In the case of certain treatments under way, with corticoids or with drugs acting on the phosphocalcic metabolism, the use of HATriC™ must be very carefully evaluated.
- When there is significant vascular impairment proximal to the graft site.
- When used with Arthrex Opening Osteotomy Systems' or substantially equivalent systems' plates and screws, refer to the instructions for use for these systems for additional contraindications.

C. ADVERSE EFFECTS

Possible adverse effects include but are not limited to:

- Infections, both deep and superficial.
- Haematoma.
- Incomplete, or lack of, osseous ingrowth into bone void, as is possible with any bone void filler.
- In the United States, when used with Arthrex Opening Osteotomy systems' or substantially equivalent systems' plates and screws, refer to the instructions for use for these systems for additional adverse effects.

D. PRECAUTIONS FOR USE

- Maximum contact between HATriC™ and recipient bone must be achieved.
- Adequate post-operative immobilization is to be used with HATriC™ to prevent movement that may lead to fibrous encapsulation.
- HATriC™ must not be used for femoral or tibial osteotomies without appropriate opening osteotomy system devices, plates and screws.
- Until healing is complete, the fixation should be considered temporary, and may not withstand weight bearing or other unsupported stresses.
- When used with Arthrex Opening Osteotomy Systems' or substantially equivalent systems' plates and screws, refer to the instructions for use for these systems for post-operative rehabilitation protocols documenting patient non-weight bearing requirements and patient immobilization requirements required for femoral or tibial opening osteotomies.
- Abrasion by a metallic surface could cause the HATriC™ to degrade.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant are important considerations in the successful utilization of this device.
- An internal fixation device must never be re-used.

E. PREPARATION

Soak the selected HATriC™ wedges in sterile water or the patient's own physiological serum.. The surgeon may also use the patient's own blood, plasma-rich protein, or bone marrow in such a way that HATriC™ is completely impregnated. Based on experiments, the ideal proportion of fluid should be between 1mL for every 2g of HATriC™. The prepared HATriC™ mixture must be used immediately to preserve cell vitality. Take care when mixing HATriC™ with soaking fluids to keep the HATriC™ whole, to not destroy the macropores' structure.

F. STERILIZATION

The HATriC™ wedges are supplied as a single use, sterile device. Sterilization is accomplished by beta irradiation. HATriC™ must not be re-sterilized by any method.

G. PACKAGING AND LABELLING

- Arthrex implants should be accepted only if the factory packaging and labelling arrive intact.
- Contact Customer Service if package has been opened or altered.

H. STORAGE AND SHELF-LIFE

The shelf life of HATriC™ product is 5 years. Products must be stored in the original unopened packaging in a dry place at ambient temperatures, and should not be used after the expiration date. Should the sterile barrier appear broken, the device is not to be used. The device is not to be used after the expiration date.

I. CONTENT

| | Small wedge | Large wedge |
|--------------|-------------|-------------|
| | | |
| AR-13000B-07 | 6mm | 7mm |
| AR-13000B-10 | 9mm | 10mm |
| AR-13000B-12 | 11mm | 12mm |
| AR-13000B-15 | 14mm | 15mm |

| | Description |
|--------------|-----------------|
| AR-13005B-03 | 2mm x 3mm wedge |
| AR-13005B-05 | 2mm x 5mm wedge |
| AR-13005B-07 | 2mm x 7mm wedge |
| AR-13005B-08 | 8mm x 8mm block |

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