<p><i>Symbols Used on Labeling</i> Auf der Verpackung verwendete Symbole Symboles utilisés sur l’étiquetage Simboli usati sull’etichetta Símbolos utilizados en las etiquetas Símbolos utilizados na embalagem</p>

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English		DFU-0149r4	
A. REFERENCES			
These instructions were developed using the guidance given in the following standards:			
<ul style="list-style-type: none"> ANSI/AAMI ST79: Comprehensive guide to steam sterilization and sterility assurance in healthcare facilities. AAMI TIR12: Designing, testing, and labeling reusable medical devices for reprocessing in healthcare facilities: A guide for device manufacturers. 			
B. CLEANING AND STERILIZATION PROCEDURE			
1. No assembly/dissassembly of these devices is required. When properly performed, cleaning and/or sterilization do not compromise the use and mechanical performance of these devices. These devices are used with or on patients who may harbor both recognized and unrecognized infections. To prevent the spread of infection, all reusable devices must be thoroughly cleaned and sterilized after use on each patient.			

C. LIMITATIONS ON REPROCESSING
Repeated processing does affect the useful life of these devices. End of useful life (which is typically 1 year from the start of use) is normally determined by wear due to repeated processing, wear due to use and / or damage.

D. VALIDATION
The following recommended cleaning and sterilization methods have been validated in compliance with federal and international guidance/standards. Cleaning and sterilizing equipment and materials vary in performance characteristics. Therefore, it is the responsibility of the facility / end user to perform the appropriate validation testing.

E. CONTAMINATION AND TRANSPORTATION
It is recommended that devices are reprocessed as soon as is reasonably practical following use.

- F. CAUTIONS**
- Never immerse the device in water or cleaning detergent.
 - Never clean the device with bleach, chlorine-based detergents, liquid or chemical disinfectants, or any products containing sodium hydroxide (such as INSTRUKLENZ or Buell Cleaner). These fluids will degrade the anodized aluminum coating and potentially damage the device's functionality.
 - Never clean the device in an ultrasonic cleaner or a combination washer/sanitizer.
 - To avoid damaging the devices, do not impact or subject the devices to blunt force. Rough handling of the devices may cause damage to the devices.
 - Do not use the device for any purpose other than its intended use. Failure to use this device in accordance with the directions for its below may result in device failure, render the device unsuitable for its intended use, or com promise the procedure.
 - Only use the APS II console that has been developed by Arthrex specifically for the handpieces.
 - Do not attach the handpiece during APS II Console Self Test or the Brightness or Language Selection Programming Modes. Refer to the *APS II User’s Guide* for more information.
 - To avoid damaging the device and disposables used with the device, do not subject the disposable to excessive force. Excessive force may render the disposable and/or device inoperable.
 - Devices with adjustable components must be handled with care. Overtightening or rough handling of the instrument may damage the locking mechanism.

- G. PREPARATION FOR CLEANING**
- Disconnect accessory cable from handpiece, where applicable.
 - Open all suction pathways and chuck openings to enable maximum exposure.
 - Remove dried-on soil from devices, especially in areas such as joints and crevices, prior to washing.

- H. INSPECTION AND MAINTENANCE**
- Arthrex devices are precision medical instruments and must be used with care.
 - Inspect the device for damage prior to use and at all stages of handling thereafter.
 - If damage is detected, consult your Arthrex representative for guidance.
 - Check device for visible soil. Clean the device using the guidelines for manual cleaning if any soil is visible.
 - Lubricate all moving parts with a soluble instrument lubricant prior to sterilization.**

- I. MANUAL CLEANING**
- Perform rinsing and cleaning as soon as possible after use.
 - Thoroughly scrub the device with a clean, soft brush that is dampened with a mild, enzymatic, pH-balanced detergent. Always avoid any harsh materials that can scratch or mar the surface of the device.
 - Remove all traces of blood, debris, and stains.
 - Do not immerse the device in soap solution or rinse-water.
- a. To clean the cannulated section of the device:
- Feed the wire end of a cleaning brush through the cannulation of the device.
 - Pull the brush completely through and repeat until all soil is removed.
6. Manipulate all moving parts of the device to ensure that all soil is removed. If any soil remains, clean the device again until the soil is removed.
7. Keep the nose of the device pointed downward and rinse it under running water to remove all traces of soap.

Note: If using tap water to clean with, flush all surfaces with distilled water to prevent metal discoloration.

- After flushing operation is complete, it is recommended that forced air be used through all channels to remove water droplets.
- Gently shake the equipment free of water and wipe the surfaces with a clean, lint-free towel.

J. STERILIZATION
Sterilization does not eliminate the need for proper cleaning of devices prior to sterilization. Sterilization capabilities, cleaning, handling, and storage of devices are the responsibility of qualified facility/user personnel.
The APS II accessory handpieces and cables may be sterilized by using one of the following autoclaving parameters:

STERILIZATION PARAMETERS: FOR THE USA ONLY:			
	Exposure Temperature	Exposure Time	Drying Time
G r a v i t y - D i s p l a c e m e n t S t e a m S t e r i l i z a t i o n C y c l e	132°C (270°F)	15 Minutes	15 to 30 Minutes
	135°C (275°F)	10 Minutes	30 Minutes
	121°C (250°F)	30 Minutes	15 to 30 Minutes
P r e - v a c u u m C y c l e	132°C (270°F)	4 Minutes	20 to 30 Minutes
	135°C (275°F)	4 Minutes	20 to 30 Minutes

CAUTION: NEVER allow the APS II console plug pins to have any contact with liquids. Remove dust or moisture, if necessary, with dry compressed air. ONLY dry connectors may be plugged into the APS II console.

Cooling – The device must be adequately cooled, after being removed from the sterilizer. It should not be touched during the cooling process. Do not place the device on a cold surface or immerse in a cold fluid.
CAUTION: NEVER submerge the cable or handpiece in liquid to cool, as this can damage the cable, handpiece, connectors, and/or seals.

- K. PACKAGING**
- A standard packaging material may be used. Ensure that the pack is large enough to contain the device.
 - Devices may be loaded into dedicated instrument trays or general-purpose sterilization trays. Ensure that the device is protected and does not exceed 8.5 kg/18.7 lb per tray. Wrap the trays using the appropriate method.

L. STORAGE
Store in a cool, dry location.

M. SPECIAL PRECAUTION - TRANSMISSIBLE SPONGIFORM ENCEPHALOPATHY AGENTS

- It is outside the scope of this document to describe in detail the precautions that should be taken for Transmissible Spongiform Encephalopathy Agents.
- The agents for transmission of Creutzfeldt-Jakob disease are believed to be resistant to normal processes of disinfection and sterilization and therefore the normal processing methods of decontamination and sterilization as described above may not be appropriate where CJD transmission is a risk.
- In general, the issues that come into contact with orthopedic surgical devices are those of low TSE infectivity. However, particular precautions should be taken when handling devices that have been used on known, suspected, or at-risk patients.

- N. DIRECTIONS FOR USE**
- Users of this device are encouraged to contact their Arthrex representative if, in their professional judgment, they require a more comprehensive surgical technique.
 - Insert the AR-8320C cable connector into the cable receptacle of the handpiece so that the red dots on the connector and receptacle align and engage easily.
 - Connect the other end of the AR-8320C cable to the handpiece receptacle of the console so that the red dots on the connector and receptacle align and engage easily.
 - The accessory cable connects and locks to the console and handpiece to prevent accidental separation during use. To avoid damage, only disconnect the accessory cable by pulling on the cable connector shell (plug).
 - The console will recognize the handpiece type, and display available options for the operator to adjust, as desired.

- I. Drill Option**
- The drill handpiece is equipped with a Jacobs chuck. A corresponding key is enclosed to attach desired accessory components, such as drill pins or reamers.
 - The speed can be adjusted in 5 steps by using the footswitch SPEED button or the RPM selection button on the front panel of the console. The corresponding setting is shown on the VFD display for the selected channel. Available speed settings are: 100, 300, 500, 900 and 1400 rpm.
 - Drill direction is selected with the FORWARD or REVERSE footswitch pedal, or with the FORWARD or REVERSE button on the front panel of the console. The OSCILLATE function is disabled for the drill.

- II. Sagittal Saw Option**
- Zimmer/Hall blades may also be used in the chuck of the sagittal saw. The blade angle can be adjusted vertically up to 180° using the interlocking mechanism and the key.
 - The sagittal saw handpiece operates at a speed of 18,000 cycles per minute. This preprogrammed setting cannot be changed.
 - Activate the saw by selecting the FORWARD or REVERSE footswitch pedal, or with the FORWARD or REVERSE button on the front panel of the console. The OSCILLATE function is disabled for the saw.

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