Symbols Used on Labeling
 Auf der Verpackung verwendete Symbole
 Symboles utilisés sur l'étiquetage
 Simboli usati sull'etichetta
 Símbolos utilizados en las etiquetas
 Símbolos utilizados na embalagem

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English DFU-0174r3

A. DEVICE DESCRIPTION

The Centerline™ Endoscopic Soft Tissue Release System is a handled device with a clear tip and a disposable blade that slides forwards and backwards. It can be attached to any standard video camera and light source.

The Guided Release System (GRS) Disposable blade is used with the Slotted Instrument Sheath, providing a single portal system. It can be attached to any standard 4 mm arthroscope.

B. INDICATIONS

The Centerline™ Endoscopic Soft Tissue Release System is indicated for patients diagnosed with carpal tunnel syndrome, Plantar fasciitis, and Equinus Contracture.

C. CONTRAINDICATIONS

1. Patients in which visualization of anatomic details fails to provide orientation to critical anatomy should not undergo endoscopic release.
2. Patients with known abnormalities of the wrist affecting the carpal tunnel, including distal radius and wrist deformities, rheumatoid, or other diseases of the synovium.
3. Patients with congenital anatomical abnormalities.
4. Patients with smaller hands or an unusually tight carpal tunnel.
5. Blood supply limitations and previous infections, which may retard healing.
6. Any active infection or blood supply limitations.
7. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
8. Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.

E. WARNINGS

1. Do not re-sterilize this device.
2. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.
3. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

F. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

G. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.

H. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

I. STORAGE CONDITIONS

This device must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

J. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing

any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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Centerline™ Endoscopic Soft Tissue Release System

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Arthrex, Inc.
Arthrex GmbH