**Symbols Used on Labeling**  
Auf der Verpackung verwendete Symbole  
Symboles utilisés sur l'etiquetage  
Simboli usati sull'etichetta  
Símbolos utilizados en las etiquetas  
Símbolos utilizados na embalagem

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**English** **DFU-0074r II**

**A. DEVICE DESCRIPTION**  
TransFix® implants are cannulated or non-cannulated pins where one end is tapered and the other end has barbs or threads. The implant with barbs is impacted into position, while the threaded implant is screwed in. The TransFix family consists of the TransFix, TransFix T3 Pin, Bio-TransFix, Bone Tendon Bio-TransFix and BioComposite TransFix implants.

**B. INDICATIONS**  
These devices are intended to be used for fixation of tissue including ligament or tendon to bone and bone tendon bone during orthopaedic procedures. Interference fixation is appropriate for surgeries of the knee. Specifically; Anterior Cruciate Ligament Repair, Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair and Iliotibial Band Tenodesis.

**C. CONTRAINDICATIONS**  
1. Insufficient quantity or quality of bone.  
2. Blood supply limitations and previous infections, which may retard healing.  
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.  
4. Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.  
5. Any active infection or blood supply limitations.  
6. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.  
7. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.  
8. Do not use for surgeries other than those indicated.

**D. ADVERSE EFFECTS**  
1. Infections, both deep and superficial.  
2. Foreign body reactions.  
3. **Bio- only:** Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.

**E. WARNINGS**  
1. An internal fixation device must never be reused.  
2. **Bio- only:** Do not re-sterilize this device.  
3. **Titanium only:** All metallic implant devices used for this surgical procedure should have the same metallurgical composition.  
4. Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand repeat bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.  
5. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.  
6. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.  
7. Detailed instructions on the use and limitations of the device should be given to the patient.  
8. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.  
9. This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating or migration in the MR environment. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.  
10. **Bio- only:** Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.  
11. For bone-tendon grafts, do not use a bone block that is smaller than 9mm in diameter.

**F. PRECAUTIONS**  
1. The BioComposite Transfix implants require driver AR-1973CD for proper insertion. Failure to use this driver could result in implant damage.  
2. Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. Or, contact your Arthrex representative for an onsite demonstration. The Arthrex website also provides detailed surgical technique information and demonstrations.

**G. PACKAGING AND LABELING**  
1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.  
2. Contact Customer Service if the package has been opened or altered.

**H. STERILIZATION**  
The device is supplied sterile. Refer to the package label for the sterilization method.  
**Titanium only:** The device can be re-sterilized following adequate cleaning if sterility is compromised, by using one of the following sterilization parameters.  
Follow your country-specific guidelines, standards, and requirements.

STERILIZATION PARAMETERS: FOR THE USA ONLY:			
Gravity-Displacement Steam Sterilization Cycle	Exposure Temperature	Exposure Time	Drying Time
	121°C (250°F)	30 Minutes	15 to 30 Minutes
	132°C (270°F)	15 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	132°C (270°F)	4 Minutes	20 to 30 Minutes
	135°C (275°F)	3 Minutes	16 Minutes

**STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:**

Gravity-Displacement Steam Sterilization Cycle	Exposure Temperature	Exposure Time	Drying Time
	132°C - 135°C (270°F - 275°F)	18 Minutes	15 to 30 Minutes
	121°C (250°F)	30 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	132°C - 135°C (270°F - 275°F)	4 Minutes	20 to 30 Minutes

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.

Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions.  
Cooling – The device must be adequately cooled after being removed from the sterilizer.

**I. MATERIAL SPECIFICATIONS**  
Refer to the package label for the materials.  
The device is made of Titanium, Poly (L-lactide), or Poly (L-Lactide-co-D, L-Lactide) and Biphasic Calcium Phosphate.

**J. STORAGE CONDITIONS**  
Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.  
Non-Bioabsorbable devices (titanium) must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

**K. INFORMATION**  
Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. Or, contact your Arthrex representative for an onsite demonstration. The Arthrex website also provides detailed surgical technique information and demonstrations.

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