# PROWICK™ SHOULDER POSTOPERATIVE DRESSING AND COLD THERAPY SYSTEM INSTRUCTIONS FOR USE (AR-I625P)

## ENGLISH

### *Description/Indications for use:*

The ProWick Shoulder Postoperative Dressing and Cold Therapy System contains an absorbent dressing for use during the first 48 hours after arthroscopic or mini-open shoulder surgery. The single size ProWick Shoulder Dressing is designed to accommodate most shoulders and is intended for absorbing body fluids that may exude from incisions after surgery. The Arthrex ProWick Dressing is indicated for the management of post-surgical wounds. The ProWick Shoulder Postoperative Dressing and Cold Therapy System contains one shoulder cover dressing, a set of five foam islands, two compression straps, one cold pack, and one set of five waterproof bandages.

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

## ENGLISH

### Instructions for use:

**1.** Unfold one Compression Strap and place it on the O.R. table such that the patient will lie on the Compression Strap on their nonsurgical side, with the strap at the level of the patient's upper chest (Figure 1). Ensure that the end with the Velcro®-like hook strips extends 20" beyond the edge of the table, in the direction the patient will face when lying in the lateral position. (If operating on a patient in a seated position, place the Compression Strap behind the patient's back with the end containing the hook strips on the non-operative side.)

**2.** Position the patient on the O.R. table so that the strap is beneath the upper chest and next to the skin (Figure 2). Before surgical skin prep, tuck the ends of the Compression Strap away from the surgical site to keep them dry.

When the operation is complete and the arthroscopic incision sites have been closed and dried per standard hospital procedure (typically suture and/or a SteriStrip™ ), the ProWick Shoulder Postoperative Dressing System can be applied.

**3.** Peel a rectangular foam island off the release liner. (If desired and per hospital procedure, topical medication to relieve postoperative pain may be placed in the center of the foam island prior to placing the island over the wound site.) Apply the foam islands parallel to the incisions, and ensure that the adhesive is flush to the skin on all sides (Figure 3).

**4.** Place the ProWick Cover Dressing on the shoulder with the label facing out so that the "wings" fully cover the foam islands placed over the most anterior and posterior sites.

Take the end of the Compression Strap in front of the patient and, with the hook strips and label facing away from patient, attach the end to the uppermost anterior wing of the ProWick Cover Dressing, covering the hook strips (Figure 4).

**5-8.** Wrap the other end of the Compression Strap (end toward patient's back) around the shoulder in a figure-8 fashion by first passing the strap from back to front over the hook strips on the

uppermost posterior wing of the dressing and down around the opposite diagonal corner of the dressing (Figure 5). Then wrap once around the base of the dressing covering the lower hook strips (Figure 6). Complete an additional wrap up the biceps and under the axilla (Figure 7). Finish by wrapping the strap back over the shoulder and attaching it to the exposed hook strips on the end of Compression Strap already secured to the anterior wing of the dressing (Figure 8). If applicable, cut off any excess strap after securing the end to the exposed hook strips.

Adjust the Compression Strap as needed so that the fit is fairly snug, as tension will release as the surgical fluids are absorbed and the swelling decreases.

The ProWick Shoulder Dressing should be removed within 48 hours after placement. See Cold Pack and Waterproof Bandage Instructions For Use. The Instructions For Use should accompany the respective products to the location where they will be applied.

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None known at this time.

### PRECAUTIONS

**CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.**

- The ProWick Shoulder Dressing should be removed within 48 hours after placement.
- Do not immerse dressing in fluid or expose to excess moisture, as dressing could become saturated.
- ProWick Shoulder Postoperative Dressing System is for SINGLE USE ONLY. DO NOT RESTERILIZE OR REUSE.

### HOW SUPPLIED

The ProWick Shoulder Postoperative Dressing and Cold Therapy System is supplied sterile and is for single use only. Do not resterilize.

### STORAGE AND HANDLING

Store in a cool, dry place and in a manner that protects the integrity of the package and sterile barrier.

### PREPARATION FOR USE

Read the Instructions For Use in their entirety prior to using the device. Inspect for damage to product or sterile barrier. Do not use if the product is damaged or the sterile barrier is compromised.

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Velcro is a registered trademark of Velcro Industries B.V.
© 2011, Arthrex Inc. All rights reserved. LA11625, rev. 4
### Arthrex, Inc.
### Arthrex GmbH