Symbols Used on Labeling

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Sterile unless the package is damaged or open. Method of sterilization - EO

Steril, solange die Verpackung ungeöffnet und unbeschädigt ist. Sterilisationsmethode - EO
Produit stérile si l'emballage n'a pas été ouvert ou endommagé. Méthode de stérilisation - EO
Il prodotto è sterile se la confezione non è aperta o danneggiata. Metodo di sterilizzazione - EO
Esteril mientras el envase no sea abierto o dañado. Método de esterilización - EO
Estéril a não ser que a embalagem esteja danificada ou aberta. Método de esterilização - EO

Sterile unless the package is damaged or open. Method of sterilization - gamma radiation

Steril, solange die Verpackung ungeöffnet und unbeschädigt ist. Sterilisationsmethode - Bestrahlung
Produit stérile si l'emballage n'a pas été ouvert ou endommagé. Méthode de stérilisation - irradiation
Il prodotto è sterile se la confezione non è aperta o danneggiata. Metodo di sterilizzazione - Raggi Gamma
Esteril mientras el envase no sea abierto o dañado. Método de esterilización - irradiación
Estéril a não ser que a embalagem esteja danificada ou aberta. Método de esterilização – radiação gama

The product meets the essential requirements of Medical Device Directive 93/42 EEC.

Das Produkt entspricht den grundlegenden Anforderungen der Richtlinie des Rates über Medizinprodukte 93/42/EWG.
Ce produit est conforme aux exigences de la directive sur les dispositifs médicaux CEE 93/42.
Il prodotto è conforme ai requisiti essenziali della Direttiva CEE 93/42 sui Dispositivi Medici.
Este producto cumple con las normas básicas de la Directiva de productos médicos. 93/42 CEE.
O produto cumpre os requisitos essenciais da Directiva de Dispositivos Médicos 93/42 EEC.

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DFU-0152r4

A. DEVICE DESCRIPTION

The device is a metal threaded or non-threaded pin, and is available in different sizes.

B. INDICATIONS

The Arthrex K-Wire is indicated for use in fixation of bone fractures, for bone reconstructions, as guide pins for insertion of other implants or implantation through the skin so that traction may be applied to the skeletal system. Examples include:

- Fixation of small bone fragments, in long bone or small bone fractures
- Arthrodesis in hand and foot surgery
- Distal or proximal metatarsal or metacarpal osteotomies
- Mono or Bi-Cortical osteotomies in the foot or hand
- Fixation of osteotomies for Hallux Valgus treatment (such as Scarf Chevron, etc.)
- Guide Wire in hip pinning procedures
- Align and reduce long bone fractures
- For use with cerclage wire/cable in treating greater trochanter fractures

C. CONTRAINDICATIONS

- Insufficient quantities or quality of bone. Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.

E. WARNINGS

- An internal fixation device must never be re-used.
- All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating or migration in the MR environment. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging.
- Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
- Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device.
- Over-flexing of the device may cause breakage or cracking.

F. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.

G. STERILIZATION

This device is provided non-sterile and can be sterilized. It must be adequately cleaned, then sterilized using one of the following sterilization parameters.

Follow your country-specific guidelines, standards, and requirements.

STERILIZATION PARAMETERS: FOR THE USA ONLY:

	Exposure Temperature	Exposure Time	Drying Time
G r a v i t y - Displacement S t e a m Sterilization Cycle	121°C (250°F)	30 Minutes	15 to 30 Minutes
	132°C (270°F)	15 Minutes	15 to 30 Minutes
	135°C (275°F)	10 Minutes	30 Minutes
P r e - v a c u u m Cycle	132°C (270°F)	4 Minutes	20 to 30 Minutes
	135°C (275°F)	3 Minutes	16 Minutes

STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:

	Exposure Temperature	Exposure Time	Drying Time
G r a v i t y - Displacement S t e a m Sterilization Cycle	132°C - 135°C (270°F - 275°F)	18 Minutes	15 to 30 Minutes
	121°C (250°F)	30 Minutes	15 to 30 Minutes
P r e - v a c u u m Cycle	132°C - 135°C (270°F - 275°F)	4 Minutes	20 to 30 Minutes

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.

Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions.

Cooling – The device must be adequately cooled, after being removed from the sterilizer.

H. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

This device is made of Stainless Steel.

I. STORAGE CONDITIONS

Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

J. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

K. DIRECTIONS FOR USE

- Insert the K-Wire into the bone.
- Seat the pin to the approximate length necessary.
- Confirm pin placement using fluoroscopy and leave it internal to the bone as in any standard metallic fixation.

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