Symbols Used on Labeling

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English DFU-0119r5

A. DEVICE DESCRIPTION

The Osteotomy Guide Assembly is a metal assembly guide designed to align guide pins for an osteotomy.

B. INTENDED USE

The Arthrex Osteotomy Guide Assembly is intended to be used in conjunction with the Arthrex Opening Wedge Osteotomy System to locate the Osteotomy site.

C. WARNINGS

- Do not drill the pin through the Arthrex Osteotomy Guide Assembly as this may result in irreparable damage to the spring mechanism of the assembly.
- Over tightening the adjustment knob may strip the knob and result in failure to maintain the desired Parallel Guide Sleeve angle.
- Do not drill the 3.0 Guide Pin out of the Osteotomy Guide Assembly as this may result in irreparable damage to the Osteotomy Guide Assembly.

D. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if package has been opened or altered.

E. STERILIZATION

This device is provided non-sterile and can be sterilized. It must be adequately cleaned, then sterilized using one of the following sterilization parameters.

Follow your country-specific guidelines, standards, and requirements.

STERILIZATION PARAMETERS: FOR THE USA ONLY:			
	Exposure Temperature	Exposure Time	Drying Time
G r a v i t y - D i s p l a c e m e n t S t e a m S t e r i l i z a t i o n C y c l e	121°C (250°F)	30 Minutes	15 to 30 Minutes
	132°C (270°F)	15 Minutes	15 to 30 Minutes
	135°C (275°F)	10 Minutes	30 Minutes
P r e - v a c u u m C y c l e	132°C (270°F)	4 Minutes	20 to 30 Minutes
	135°C (275°F)	3 Minutes	16 Minutes

STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:			
	Exposure Temperature	Exposure Time	Drying Time
G r a v i t y - D i s p l a c e m e n t S t e a m S t e r i l i z a t i o n C y c l e	132°C - 135°C (270°F - 275°F)	18 Minutes	15 to 30 Minutes
	121°C (250°F)	30 Minutes	15 to 30 Minutes
	132°C - 135°C (270°F - 275°F)	4 Minutes	20 to 30 Minutes

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities”, for specific information.

Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer’s instructions.

Cooling – The device must be adequately cooled after being removed from the sterilizer. Do not touch the device during the cooling process. Do not place the device on a cold surface or immerse it in a cold fluid.

F. MATERIAL SPECIFICATIONS

The device is made from stainless steel and aluminum.

G. STORAGE CONDITIONS

This device must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

H. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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