## DFU-0042 Revision 7

CE  
0086

Arthrex®

### IMPORTANT PRODUCT INFORMATION

WICHTIGE PRODUKTINFORMATION

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English	DFU-0042r7
<b>A. DEVICE DESCRIPTION</b> The Ankle Distraction Strap is an adjustable strap used for effective traction and grip. The Ankle Arthroscopy Distractor Strap is used in conjunction with the Ankle Distractor Set. The adjustable strap has soft nonslip foam pads and provides effective traction and grip.	<b>K. INFORMATION</b> Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

<b>B. INDICATIONS</b> The Ankle Straps are intended to provide safe, non-invasive distraction and rotational control of the foot for ankle arthroscopy.	<b>C. CONTRAINDICATIONS</b> 1. Blood supply limitations and previous infections, which may retard healing. 2. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
--	--

<b>D. ADVERSE EFFECTS</b> 1. Allergies and other reactions to device materials. 2. Foreign body reactions.
--

<b>E. WARNINGS</b> 1. Do not re-sterilize this device. 2. Do not re-use this device. 3. Do not disassemble prior to use. 4. Do not over tighten straps. 5. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.
---

<b>F. PRECAUTIONS</b> Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
--

<b>G. PACKAGING AND LABELING</b> 1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact. 2. Contact Customer Service if the package has been opened or altered.
---

<b>H. STERILIZATION</b> The device is provided sterile. Refer to the package label for the sterilization method.
---

<b>I. MATERIAL SPECIFICATIONS</b> <b>AR-1710:</b> Cotton straps, with chrome-plated steel and nylon rings. It has an adjustable Velcro® strap. The spring clip is zinc plated steel. <b>OR</b> <b>AR-1712:</b> Nylon straps and Neoprene foam with a Silastic® medical adhesive silicone. The universal hook is made of plastic.
---

<b>J. STORAGE CONDITIONS</b> This device must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.
---

<b>Symbols Used on Labeling</b> Auf der Verpackung verwendete Symbole Symboles utilisés sur l'étiquetage Simboli usati sull'etichetta Símbolos utilizados en las etiquetas Símbolos utilizados na embalagem
--

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