Symbols Used on Labeling
Auf der Verpackung verwendete Symbole
Symboles utilisés sur l'étiquetage
Simboli usati sull'etichetta
Símbolos utilizados en las etiquetas
Símbolos utilizados na embalagem

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English

DFU-0166r6

A. INDICATIONS
The *Arthrex SportKnee™* is indicated for partial replacement of the articulating surfaces of the knee when the medial compartment is affected due to compartmental primary degenerative or post-traumatic degenerative disease. The femoral component is intended to be used with or without bone cement. The tibial component is intended to be used with bone cement.

B. CONTRAINDICATIONS

1. Symptomatic degenerative disease of the ipsilateral posterior femoral condyle, contralateral compartment or patellofemoral compartment of the same knee.
2. Instability caused by ligament deficiency.
3. Deformity or malalignment of the knee.
4. Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
5. Blood supply limitations and previous infections, which may tend to retard healing.
6. Any active infection or blood supply limitations.
7. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
8. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb, or disrupt the growth plate.
9. Pathological conditions, such as insufficient quantity or quality of bone (e.g., cystic changes or severe osteopenia), which may compromise implant fixation.
10. Do not use for surgeries other than those indicated.
11. Do not use for revision surgeries.

C. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Foreign body reactions.
3. Loosening, migration, or loss of fixation of the implant.
4. Wear and damage to the implant articulating surface.
5. Particulate wear debris and mild tissue discoloration from metallic components have been noted with other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
6. Wear and damage to the adjacent and opposed articular cartilage surfaces or soft tissue support structures.
7. Postoperative pain or incomplete resolution of preoperative symptoms.
8. Periarticular calcification or ossification with or without impeding of joint mobility.
9. Incomplete range of motion due to improper selection or positioning of components.
10. Fatigue fracture of the implants.
11. Intraoperative or postoperative bone fracture.
12. Allergic-like reactions to materials.
13. Compartment syndrome.
14. Possible necessity for revision.

D. WARNINGS

1. Do not re-sterilize the ultra high molecular weight polyethylene tibial component.
2. All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
3. Proper patient selection is critical to outcomes. Use of the SportKnee on patients with contraindications may impair implant function or lead to poor outcomes.
4. The postoperative regimen prescribed by the physician should be strictly followed.
5. Detailed instructions on the use and limitations of the device should be given to the patient.
6. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
7. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device. Patient sensitivity to the device materials should be considered prior to implantation.
8. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
9. The SportKnee has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. The SportKnee has not been tested for heating or migration in the MR environment. *The femoral component of the SportKnee is manufactured from a metallic material, and surgeons can expect that MR artifacts will be present during routine MR Imaging.*

E. PRECAUTIONS
Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print format, video format, and electronically. Or, contact your Arthrex representative for an onsite demonstration. The Arthrex website also provides detailed surgical technique information and demonstrations.

Tibial Bone Preparation

1. Extend the knee to approximately 45-60 degrees to avoid damaging the posterior condyle.
2. In order to assist with perpendicular placement of the Tibial Guide Pin in relation to the tibial plateau, verify guide placement intraarticularly and extraarticularly. Extraarticularly, the alignment rod will be parallel with the tibia in the sagittal plane when the Drill Sleeve is perpendicular to the tibial plateau.
3. To ensure accurate placement of the Tibial Guide Pin, prevent skiving of the Tibial Reamer by focusing on the outer cortex until the reamer is in the bone about 5-8 mm.
4. Use a Parallel Drill Sleeve to correct placement of the pin.

Femoral Bone Preparation

1. The shaft of the femoral drill guide must be parallel with the long axis of the femur on the sagittal plane.

Femoral Trialing

1. The Femoral Trial must be flush to 0.5 mm recessed.

2. Extreme care is needed when reaming without a positive stop.

Completion of Tibial Bone Preparation

1. Avoid contact of the tibial Retroreamer blade with the femoral condyle and tibial spine. *Femoral condyle contact:* Rotate the tibial RetroReamer by hand 360 degrees and reposition the knee. *Tibial spine contact:* Flatten the tibial spine with a burr.
2. Do not operate the tibial RetroReamer in reverse.
3. Do not over-ream the tibial socket. Use gentle pressure and check depth often with the tibial depth probe.

Implantation

1. Remove any cement that extrudes into the joint space.

F. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.

G. STERILIZATION
The device is provided sterile. Refer to the package label for the sterilization method.

If the sterility of the femoral component is compromised, first verify that the articular surface is free of damage and then use the sterilization parameters below.

NOTE: UHMWPE components may be damaged if resterilized.

Follow your country-specific guidelines, standards, and requirements.

STERILIZATION PARAMETERS: FOR THE USA ONLY:

Gravity-Displacement	Exposure Temperature	Exposure Time	Drying Time
Steam	121°C (250°F)	60 Minutes	15 to 30 Minutes
Sterilization Cycle	132°C (270°F)	5 Minutes	20 to 30 Minutes
Pre-vacuum Cycle	135°C (275°F)	5 Minutes	20 to 30 Minutes

STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:

Gravity-Displacement	Exposure Temperature	Exposure Time	Drying Time
Steam	132°C – 135°C (270°F – 275°F)	18 Minutes	15 to 30 Minutes
Sterilization Cycle	121°C (250°F)	60 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	132°C - 135°C (270°F - 275°F)	5 Minutes	20 to 30 Minutes

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST46, "Good Hospital Practice: Steam Sterilization and Sterility Assurance," for specific information.

Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions.

Cooling – The device must be adequately cooled after being removed from the sterilizer.

H. MATERIAL SPECIFICATIONS
Refer to the package label for the materials.

Femoral Implant: Cast Cobalt Chrome
Femoral Fixation Surface Coating: Plasma spray w/CP Titanium powder
Tibial Implant: Ultra High Molecular Weight Polyethylene (UHMWPE)
Tibial Suture Pin: Titanium

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**Arthrex, Inc.**
Arthrex GmbH