<p><b><i>Symbols Used on Labeling</i></b>  Auf der Verpackung verwendete Symbole  Symboles utilisés sur l'étiquetage  Simboli usati sull'etichetta  Símbolos utilizados en las etiquetas  Símbolos utilizados na embalagem</p>
--

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English	DFU-0165r4
<b>A. DEVICE DESCRIPTION</b> The Double Syringe Autologous Conditioned Plasma (ACP®) System is comprised of an inner syringe inside an outer syringe. The inner syringe was specifically designed to fit within and function in concert with the specifically designed outer syringe. The two syringes are connected via a luer lock connector. A threaded screw cap connected to the luer tip of the outer syringe can be removed prior to use. The Double Syringe Autologous Conditioned Plasma (ACP) System was designed to be used with the existing Rotofix 32 A centrifuge.	Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
<b>B. INDICATIONS</b> The Double Syringe Autologous Conditioned Plasma (ACP) System is used to facilitate the safe and rapid preparation of autologous platelet-rich plasma (PRP) from a small sample of blood at the patient's point of care. The PRP can be mixed with autograft and allograft bone prior to application to an orthopedic surgical site as deemed necessary by the clinical use requirements.	The PRP can be mixed with autograft and allograft bone prior to application to an orthopedic surgical site as deemed necessary by the clinical use requirements.
<b>C. CONTRAINDICATIONS</b> 1. Blood supply limitations and previous infections, which may retard healing. 2. Any active infection or blood supply limitations.	
<b>D. ADVERSE EFFECTS</b> 1. Infections, both deep and superficial. 2. Allergies and other reactions to device materials. 3. Hematoma. 4. Damage to blood vessels and nerve damage resulting in pain or numbness from autologous sampling. 5. Delayed wound healing.	
<b>E. WARNINGS</b> 1. Do not re-sterilize this device. 2. Do not re-use this device. 3. Use only the Rotofix 32 A Benchtop Centrifuge, and follow the manufacturer's instructions. Outcomes using centrifuges from other manufacturers are unknown. 4. Follow the manufacturer's directions that are provided with the package insert for the Anticoagulant Citrate Dextrose Solution A (ACD-A). 5. The safety and effectiveness of this device for bone healing and hemostasis have not been established. 6. The PRP must be used within four (4) hours of blood collection. 7. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device. 8. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.	
<b>F. PACKAGING AND LABELING</b> 1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact. 2. Contact Customer Service if the package has been opened or altered.	
<b>G. STERILIZATION</b> The device is supplied sterile. Sterilization is accomplished by Ethylene Oxide sterilization. Serious allergic reactions to Ethylene Oxide sterilization have been noted in some sensitized individuals. Certain Arthrex devices that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.	
<b>H. MATERIAL SPECIFICATIONS</b> Refer to the package label for the materials. <b>Outer Syringe:</b> Polypropylene, TPV, Silicone Oil <b>Inner Syringe:</b> Polypropylene, Silicone Oil <b>Cap:</b> Polyethylene	
<b>I. STORAGE CONDITIONS</b> These devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.	
<b>J. INFORMATION</b> Surgeons are advised to review the product specific surgical technique prior to performing any surgery.	

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**Arthrex, Inc.**
**Arthrex GmbH**