Symbols Used on Labeling

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Símbolos utilizados na embalagem

All symbols shown may not apply to this device. Please refer to the package labeling for utilized symbols.

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

English DFU-0039r12

A. DEVICE DESCRIPTION
TissueButtons are bioabsorbable buttons that are available in several sizes. Each button has four holes that can be used for suture fixation.

B. INDICATIONS
The TissueButton™ is intended for suture fixation of soft tissue to bone, utilizing the Arthrex Corkscrew® Suture Anchor and polyester braided sutures appropriate to the procedure, up to size 5 (single strand). This product is intended for use in repairing rotator cuff repairs.

- C. CONTRAINDICATIONS**
1. Insufficient quantity or quality of bone.
 2. Blood supply limitations and previous infections that may tend to retard healing.
 3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
 4. Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
 5. Any active infection or blood supply limitations.
 6. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
 7. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
 8. Do not use for surgeries other than those indicated.

- D. ADVERSE EFFECTS**
1. Infection, both deep and superficial.
 2. Foreign body reactions.
 3. Allergic –like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.

- E. WARNINGS**
1. An internal fixation device must never be re-used.
 2. Do not resterilize this device.
 3. Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
 4. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
 5. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
 6. Detailed instructions on the use and limitations of this device should be given to the patient.
 7. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
 8. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
 9. Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.

- F. PACKAGING AND LABELING**
1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
 2. Contact Customer Service if package has been opened or altered.

G. STERILIZATION
This device is provided sterile. Refer to the package label for the sterilization method.
Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.

H. MATERIAL SPECIFICATION
Refer to the package label for the materials.
This device is made of Poly (L-Lactide).

I. STORAGE CONDITIONS
Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.

J. INFORMATION
Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical

technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

K. DIRECTIONS FOR USE
Users of this device are encouraged to contact their Arthrex representatives if, in their professional judgment, they require a more comprehensive surgical technique.

1. Suture is passed through the soft tissue to be repaired according to surgeon's preference.
2. The two free ends of the suture are brought through the TissueButton in two different holes.
3. Soft tissue reattachment is accomplished by tying the suture ends over the TissueButton, securing the suture and button to the bone.

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TissueButton™ DFU-0039

Revision 12

TissueButton
TissueButton
TissueButton
TissueButton
TissueButton

IMPORTANT PRODUCT INFORMATION
WICHTIGE PRODUKTINFORMATION
NOTICE D'UTILISATION IMPORTANTE
IMPORTANTI INFORMAZIONI PER L'USO
INSTRUCCIONES IMPORTANTES PARA EL USO
INFORMAÇÕES IMPORTANTES ACERCA DO PRODUTO

 **Arthrex, Inc.**
1370 Creekside Blvd.
Naples, FL 34108-1945 • USA
Toll free: 1-(800) 934-4404
www.arthrex.com

Arthrex GmbH
Erwin-Hielscher-Strasse 9
81249 München, Germany
Tel: +49 89 909005-0
www.arthrex.de

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