<p>Symbols Used on Labeling Auf der Verpackung verwendete Symbole Symboles utilisés sur l'etiquetage Simboli usati sull'etichetta Símbolos utilizados en las etiquetas Símbolos utilizados na embalagem</p>

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English	DFU-0104r6
A. DEVICE DESCRIPTION The Viper products are hand held instruments used for passing suture in an open or arthroscopic procedure.	the bottom jaw. Push the trigger to the forward position so that the suture is retracted back into the lower jaw. It is recommended that the #2 suture strand is loaded approximately 2 cm from the end of the suture strand.
B. INDICATIONS The Arthrex Viper Suture Passer and Bankart Viper Suture Passer are intended to pierce tissue with a hooked needle and retrieve a loop of suture from the other side.	<ol style="list-style-type: none"> Insert the closed jaws of the Viper through an appropriately sized cannula. Before opening the instrument jaws, ensure that the entire jaw has cleared the distal end of the cannula. The cannula can be backed up to provide more room for this maneuver. The Arthrex Notched Cannula is recommended for the larger Viper Suture Passer. Grasp tissue with the jaws of the instrument. It is recommended that 2 mm of free space remains between the edge of the tissue and the crutch of the instrument jaws. Pierce the tissue and firmly hold the handles completely closed. Whenever possible, rotate the tip of the instrument and position the scope to visualize the bottom of the lower instrument jaw to confirm that the tip of the needle has passed through the tissue and is visible in the small window. Pull the trigger to the back position to transfer the suture into the hook of the needle. Open the handles while pushing the instrument medial to the tissue edge. In cases where the tissue is highly mobile, an accessory grasping instrument or traction stitch may be necessary to control the tissue. The built-in spring arm will help hold the tissue down while pulling the suture through. After the needle and suture have been pulled through the tissue, pull the instrument back beyond the edge of the tissue and close the handles. It may be necessary to reposition the cannula to be able to clear the edge of the tissue before closing the jaws. Pull the instrument, with the loop of suture, out of the cannula.
C. ADVERSE EFFECTS <ol style="list-style-type: none"> Infections, both deep and superficial. Foreign body reactions. 	
D. WARNINGS <ol style="list-style-type: none"> The Arthrex Viper Suture Passers are precision medical instruments and must be used with care. Prior to use, inspect and verify that the device works properly. The needle must be completely through the tissue and held there while transferring the suture. The needle may not be able to penetrate completely through extremely thick tissue. CAUTION: Over stressing the handles in an attempt to fully penetrate extremely thick tissue can result in instrument breakage and loss of function. Do not pull back on the instrument while opening the handles. Carefully push forward while opening. Do not force the instrument if it seems to be stuck. Slightly wiggle and push forward while opening to clear. CAUTION: Wrenching the instrument around in an attempt to dislodge it can result in the needle bending or breaking. If the needle cannot be gently freed from the tissue, use a grasping instrument through another portal to help free the needle. Do not use to manipulate tissue or for any use other than indicated. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure. 	
E. PACKAGING AND LABELING <ol style="list-style-type: none"> Arthrex devices should be accepted only if the factory packaging and labeling arrive intact. Contact Customer Service if the package has been opened or altered. 	
F. STERILIZATION This device is provided non-sterile and can be sterilized. It must be adequately cleaned, then sterilized using one of the following sterilization parameters. <ol style="list-style-type: none"> Follow your country-specific guidelines, standards, and requirements. 	

STERILIZATION PARAMETERS: FOR THE USA ONLY			
	Exposure Temperature	Exposure Time	Drying Time
Gravity-Displacement Steam Sterilization Cycle	121°C (250°F)	30 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	132°C (270°F)	4 Minutes	20 to 30 Minutes
	135°C (275°F)	3 Minutes	16 Minutes

STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:			
	Exposure Temperature	Exposure Time	Drying Time
Gravity-Displacement Steam Sterilization Cycle	132°C – 135°C (270°F – 275°F)	18 Minutes	15 to 30 Minutes
Pre-vacuum Cycle	121°C (250°F)	30 Minutes	15 to 30 Minutes
	132°C - 135°C (270°F - 275°F)	4 Minutes	20 to 30 Minutes

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information. Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions. Cooling – The device must be adequately cooled, after being removed from the sterilizer.

G. MATERIAL SPECIFICATIONS			
This device is made out of stainless steel. Refer to the package label for the materials.			
H. STORAGE CONDITIONS			
Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.			

I. INFORMATION			
Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.			
J. DIRECTIONS FOR USE			
Users of this device are encouraged to contact their Arthrex representatives if, in their professional judgment, they require a more comprehensive surgical technique.			
WARNING: IT IS THE SURGEON'S RESPONSIBILITY TO BECOME FAMILIAR WITH THE USE AND FUNCTION OF THE VIPER BEFORE IT IS USED ON A PATIENT.			
<ol style="list-style-type: none"> To load a #2 suture: Open the instrument jaws and pull the trigger to the back position. Place the suture in the open slot on 			

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Arthrex, Inc.
Arthrex GmbH