<p>Symbols Used on Labeling</p> <p>Auf der Verpackung verwendete Symbole</p> <p>Simboli utilizzati sur l'etiquetage</p> <p>Simboli usati sull'etichetta</p> <p>Símbolos utilizados en las etiquetas</p> <p>Símbolos utilizados na embalagem</p>
--

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English	DFU-0173r5
A. DEVICE DESCRIPTION The Arthrex GraftBolt® device is a cannulated screw, with a mating sheath. These devices are available in different sizes. The sheath has ribbed tabs that deploy as the mating screw advances within the sheath.	2. Contact Customer Service if the package has been opened or altered.
B. INDICATIONS The GraftBolt is intended to be used for fixation of tissue including ligament or tendon to bone and bone tendon bone during cruciate ligament reconstruction procedures.	H. STERILIZATION This device is provided sterile. Refer to the package label for the sterilization method. Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.
C. CONTRAINDICATIONS 1. Insufficient quantity or quality of bone. 2. Blood supply limitations and previous infections, which may retard healing. 3. Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation. 4. Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions. 5. Any active infection or blood supply limitations. 6. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period. 7. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate. 8. Do not use for surgeries other than those indicated.	I. MATERIAL SPECIFICATIONS Refer to the package label for the materials. The screw and sheath are made of Polyetheretherketone (PEEK). The Bio-composite screw and sheath are made of Poly (L, DL-lactide; PLDLA), biphasic tricalcium phosphate (TCP) and hydroxyapatite (HA).
D. ADVERSE EFFECTS 1. Infections, both deep and superficial. 2. Foreign body reactions. 3. Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.	J. STORAGE CONDITIONS Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date. Non-Bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.
E. WARNINGS 1. The GraftBolt implant consists of a prepackaged mating sheath and screw to avoid mis-mating sizes. If a mating sheath or screw is compromised, then a new GraftBolt implant (sheath and screw) of the same size must be opened to obtain a new sheath or screw. Mismatched GraftBolt sheath and screw can jeopardize the fixation. 2. No other Arthrex or competitor sheath or screw is to be used with either the GraftBolt sheath or screw. Mismatched implants can jeopardize the fixation. 3. An internal fixation device must never be re-used. 4. Do not re-sterilize this device. 5. Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device. 6. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device. 7. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. To remove the GraftBolt Arthrex recommends the use of the GraftBolt Removal Tool. Device removal should be followed by adequate postoperative management. 8. Detailed instructions on the use and limitations of the device should be given to the patient. 9. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user. 10. This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating or migration in the MR environment. If the implant is manufactured from a metallic material, surgeons can expect that MR artifacts will be present during routine MR imaging. 11. Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Migration of implant position resulting in injury; (2) Risk of additional injury from postoperative trauma; (3) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (4) Pain, discomfort, or abnormal sensations due to the presence of the device; (5) Possible increased risk of infection; and (6) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture. 12. Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.	K. INFORMATION Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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**Arthrex, Inc.**
Arthrex GmbH