<p><b>Symbols Used on Labeling</b></p> <p>Auf der Verpackung verwendete Symbole</p> <p>Symboles utilisés sur l'etiquetage</p> <p>Simboli usati sull'etichetta</p> <p>Símbolos utilizados en las etiquetas</p> <p>Símbolos utilizados na embalagem</p>
---

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<b>English</b>	<b>DFU-0095r8</b>																					
<p><b>A. DEVICE DESCRIPTION</b> The NeedlePunch® device is a handled instrument. An Arthrex needle with suture is loaded into the bottom jaw of the device and it's deployed by squeezing the handle. Arthrex needles for the NeedlePunch are available in several different configurations.</p> <p><b>B. INDICATIONS</b> Arthrex NeedlePunch devices are intended to grasp tissue and pass a specific needle, with suture attached, through the tissue where it is captured and retrieved.</p> <p><b>C. WARNINGS</b></p> <ol style="list-style-type: none"> <li>Use only Arthrex needles specially designed for use with these devices.</li> <li>Prior to use, inspect and verify that the device works properly.</li> <li>The needle must be properly loaded in the lower jaw before it can be advanced through tissue with the Pushrod.</li> <li>Do not use with a damaged or bent needle. Do not force the needle into the slot of the NeedlePunch.</li> <li>Push directly in line with the Pushrod shaft when advancing the needle. Do not bend the Pushrod.</li> <li>Do not twist the punch while removing. This may make the needle disengage from the up- per jaw.</li> <li>If the needle disengages from the upper jaw capture mechanism, pull the needle forward through the tissue with a grasping instrument.</li> <li>Do not pull the needle backward by the suture. This could cause the needle to separate from the suture and result in a lost needle.</li> <li>Remove the Pushrod for cleaning and sterilization.</li> <li>Arthrex.NeedlePunch devices should be properly cleaned and sterilized to minimize the health and/or safety risks of cross-infection.</li> <li>Do not use for surgeries other than those indicated</li> <li>Do not follow these guidelines may result in instrument failure.</li> </ol>																						
<p><b>D. PRECAUTIONS</b></p> <ol style="list-style-type: none"> <li>Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.</li> </ol>																						
<p><b>E. PACKAGING AND LABELING</b></p> <ol style="list-style-type: none"> <li>Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.</li> <li>Contact Customer Service if the package has been opened or altered.</li> </ol>																						
<p><b>F. STERILIZATION</b> This device is provided sterile or non-sterile. Refer to the package label for the sterilization method. <b>Needles and Suture can't be re-sterilized.</b> <b>Devices made of metal only:</b> This device can be resterilized. It must be adequately cleaned, then sterilized using one of the following sterilization parameters. Follow your country-specific guidelines, standards, and requirements.</p>																						
<p><b>STERILIZATION PARAMETERS: FOR THE USA ONLY:</b></p> <table border="1"> <thead> <tr> <th></th> <th>Exposure Temperature</th> <th>Exposure Time</th> <th>Drying Time</th> </tr> </thead> <tbody> <tr> <td rowspan="3"><b>Gravity-Displacement Steam Sterilization Cycle</b></td> <td>121°C (250°F)</td> <td>30 Minutes</td> <td>15 to 30 Minutes</td> </tr> <tr> <td>132°C (270°F)</td> <td>15 Minutes</td> <td>15 to 30 Minutes</td> </tr> <tr> <td>135°C (275°F)</td> <td>10 Minutes</td> <td>30 Minutes</td> </tr> <tr> <td rowspan="2"><b>Pre-vacuum Cycle</b></td> <td>132°C (270°F)</td> <td>4 Minutes</td> <td>20 to 30 Minutes</td> </tr> <tr> <td>135°C (275°F)</td> <td>3 Minutes</td> <td>16 Minutes</td> </tr> </tbody> </table>			Exposure Temperature	Exposure Time	Drying Time	<b>Gravity-Displacement Steam Sterilization Cycle</b>	121°C (250°F)	30 Minutes	15 to 30 Minutes	132°C (270°F)	15 Minutes	15 to 30 Minutes	135°C (275°F)	10 Minutes	30 Minutes	<b>Pre-vacuum Cycle</b>	132°C (270°F)	4 Minutes	20 to 30 Minutes	135°C (275°F)	3 Minutes	16 Minutes
	Exposure Temperature	Exposure Time	Drying Time																			
<b>Gravity-Displacement Steam Sterilization Cycle</b>	121°C (250°F)	30 Minutes	15 to 30 Minutes																			
	132°C (270°F)	15 Minutes	15 to 30 Minutes																			
	135°C (275°F)	10 Minutes	30 Minutes																			
<b>Pre-vacuum Cycle</b>	132°C (270°F)	4 Minutes	20 to 30 Minutes																			
	135°C (275°F)	3 Minutes	16 Minutes																			
<p><b>STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY:</b></p> <table border="1"> <thead> <tr> <th></th> <th>Exposure Temperature</th> <th>Exposure Time</th> <th>Drying Time</th> </tr> </thead> <tbody> <tr> <td rowspan="2"><b>Gravity-Displacement Steam Sterilization Cycle</b></td> <td>132°C - 135°C (270°F - 275°F)</td> <td>18 Minutes</td> <td>15 to 30 Minutes</td> </tr> <tr> <td>121°C (250°F)</td> <td>30 Minutes</td> <td>15 to 30 Minutes</td> </tr> <tr> <td><b>Pre-vacuum Cycle</b></td> <td>132°C - 135°C (270°F - 275°F)</td> <td>4 Minutes</td> <td>20 to 30 Minutes</td> </tr> </tbody> </table>			Exposure Temperature	Exposure Time	Drying Time	<b>Gravity-Displacement Steam Sterilization Cycle</b>	132°C - 135°C (270°F - 275°F)	18 Minutes	15 to 30 Minutes	121°C (250°F)	30 Minutes	15 to 30 Minutes	<b>Pre-vacuum Cycle</b>	132°C - 135°C (270°F - 275°F)	4 Minutes	20 to 30 Minutes						
	Exposure Temperature	Exposure Time	Drying Time																			
<b>Gravity-Displacement Steam Sterilization Cycle</b>	132°C - 135°C (270°F - 275°F)	18 Minutes	15 to 30 Minutes																			
	121°C (250°F)	30 Minutes	15 to 30 Minutes																			
<b>Pre-vacuum Cycle</b>	132°C - 135°C (270°F - 275°F)	4 Minutes	20 to 30 Minutes																			
<p>Certain Arthrex devices that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information. Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions. Cooling – The device must be adequately cooled after being removed from the sterilizer. Do not touch the device during the cooling process. Do not place the device on a cold surface or immerse it in a cold fluid.</p>																						
<p><b>G. MATERIAL SPECIFICATIONS</b> Refer to the package label for the materials. NeedlePunch devices and needles are made of stainless steel. FiberWire® suture is made of Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester braided over a UHMWPE core. Some suture may be coreless. The UHMWPE suture is braided from 100% UHMWPE. Additional materials to the Arthrex Suture family include silicone elastomer coating, and cyanoacrylate. Suture strands that are dyed black are made of nylon. The coating acts as a lubricant for suture sliding, knot tying, and ease of passing suture through tissue. The Arthrex Suture family is available in several sizes. FiberWire® sutures meet U.S. and European Pharmacopeia standards for non-absorbable suture, except diameter, and UHMWPE sutures meet USP standards. The Arthrex Suture is supplied sterile in pre-cut lengths, in loop configurations, with swaged needles. The Arthrex Suture is available with stiffened ends. The Arthrex Suture family is available</p>																						

non-dyed, dyed, and fully or partially striped. Dyes may include: D&C Blue No. 6, D&C Green No. 6, and Logwood Black.

**H. STORAGE CONDITIONS**  
Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date. Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

**I. INFORMATION**  
Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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