Symbols Used on Labeling
Auf der Verpackung verwendete Symbole
Symboles utilisés sur l'étiquette
Simboli usati sull'etichetta
Símbolos utilizados en las etiquetas
Símbolos utilizados na rotulagem

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English

A. DEVICE DESCRIPTION
The iBalance® PFJ consists of femoral components and patella components. All components are available in a range of sizes to fit varying anatomical requirements. Femoral components are available in left and right varieties. Patella components have a circular profile. Any iBalance PFJ patella component can be used with any iBalance PFJ femoral component.

B. INDICATIONS
Indicated for use in patellofemoral knee arthroplasty in patients with:
1. Degenerative arthritis in the trochlea and patella.
2. A history of patella dislocation or fracture.
3. Failed previous surgery (arthroscopy, tibial tubercle elevation, lateral release) where pain, deformity or dysfunction persists.

These components are single use only and are intended for implantation with bone cement.

When used concurrently, the Arthrex iBalance UKA and PFJ systems create the **Arthrex iBalance BiCompartmental Arthroplasty System**. The **Arthrex iBalance BiCompartmental Arthroplasty System** is intended to be used as a multi-compartmental knee arthroplasty in patients with:

- Moderately disabling joint disease of the knee resulting from painful osteoarthritis or post traumatic arthritis;
- Correction of functional deformities;
- Revision of previous unsuccessful partial knee replacement or other procedure.

The BiCompartmental Arthroplasty System is not intended to be used as a dual-condyle or tri-compartmental knee.

These components are single use only and are intended for implantation with bone cement.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- Do not use for surgeries other than those indicated.
- An overweight or obese patient can produce loads on the prosthesis which can lead to failure of fixation or failure of the device itself.
- Severe deformity and/or recurrent subluxation of the knee joint.
- Inflammatory arthritis conditions, particularly rheumatoid arthritis.
- Ligamentous instability or soft-tissue laxity such that the postoperative stability afforded by the prosthesis may be compromised or dislocate.

D. ADVERSE EFFECTS

- Loosening, migration, or loss of fixation of the implant.
- Foreign body reactions.
- Infections, both deep and superficial.
- Particulate wear debris and mild tissue discoloration from metallic components have been noted with other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
- Wear and damage to the articulating surfaces or soft tissue support structures.
- Postoperative symptoms include, but are not limited to: pain, persistent swelling, stiffness, limited range of motion, or incomplete resolution of preoperative symptoms.
- Periarticular calcification or ossification with or without impediment of joint mobility.
- Incomplete range of motion due to improper selection or positioning of components.
- Fatigue fracture of the implants.
- Intraoperative or postoperative bone fracture.
- Allergic reactions to materials.
- Other potential adverse effects of knee surgery may include: locking, neurovascular damage, dislocation, bending of the components, leg length discrepancies, and other less common adverse effects.
- Serious complications may be associated with any total joint replacement surgery. These complications include, but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Particles from metal and polyethylene components can be shed from the components during normal use and over time. These particles may be incorporated into the joint through natural biological processes or migrate to other parts of the body. The long-term effects of these particles may include, but are not limited to, cancer, lymphadenopathy and accumulation in other tissues and organs, and systemic disease.

E. WARNINGS

- Care should be utilized in the handling of the components to minimize contamination and damage of the component surfaces.
- When using cement for fixation, the surgeon should use care to ensure complete cement contact and coverage on all parts of the prosthesis embedded in bone cement. Remove any excess cement to ensure that no cement remains on the articulating surfaces of the prosthesis.
- An internal fixation device must never be re-used.
- The surgeon must not allow damage to polished bearing surfaces because this may accelerate wear of the components. Any alteration or damage to a component may result in failure under load. Any prostheses so damaged must not be used.
- Components of the iBalance PFJ System should not be used with those of another manufacturer since articular and dimensional compatibility cannot be assured.
- Detailed instructions on the use and limitations of this device should be given to the patient. Patient conditions can affect results and outcomes. Performing activities that increase stresses on implants such as running, lifting, skiing, etc. can result in early failure of these implants.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Postoperatively and until healing is complete, the fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.

- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- This product's useful life will vary between patients, and is dependent upon accumulated weight-bearing activities and lifestyle.
- Physicians should carefully assess the activity level of the patient prior to performing arthroplasty of the knee. Increased activity as well as increased weight can lead to accelerated wear of UHMWPE components.

F. MRI SAFETY INFORMATION

1. MR Conditional
Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex iBalance PFJ is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only.
- Maximum spatial gradient magnetic field of 2000 Gauss/cm or less.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system.
- Under the scan conditions defined, the Arthrex iBalance PFJ is expected to produce a maximum temperature rise of 2 °C after 15-minutes of continuous scanning.
- Patient knee must be in direct contact with the patient table. Pillows should not be placed under the knee during scanning.

2. Artifact Information
In non-clinical testing, the image artifact caused by the Arthrex iBalance PFJ extends approximately 20 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Proper selection of the Arthrex iBalance PFJ System depends on the judgment of the surgeon.
- For cemented use only, unless otherwise indicated.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

H. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.

2. Contact Customer Service if the package has been opened or altered.

I. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

NOTE: Ultra High Molecular Weight Polyethylene (UHMWPE) components may be damaged if resterilized.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023 and ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities" for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

Femoral components are made of Cobalt Chrome.
Patella components are made of Ultra High Molecular Weight Polyethylene (UHMWPE).

K. STORAGE CONDITIONS

This device must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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iBalance® PFJ

iBalance PFJ
iBalance PFJ
iBalance PFJ
iBalance PFJ
iBalance PFJ

DFU-0183-2

Rev. 0 01/2018

IMPORTANT PRODUCT INFORMATION
WICHTIGE PRODUKTINFORMATION
NOTICE D'UTILISATION IMPORTANTE
IMPORTANTI INFORMAZIONI PER L'USO
INSTRUCCIONES IMPORTANTES PARA EL USO
INFORMAÇÕES IMPORTANTES SOBRE O PRODUTO

 **Arthrex, Inc.**
1370 Creekside Blvd.
Naples, FL 34108-1945 • USA
Toll free: 1-(800) 934-4404
www.arthrex.com

 **Arthrex GmbH**
Erwin-Hielscher-Strasse 9
81249 München, Germany
Tel: +49 89 909005-0
www.arthrex.de

This is not a warranty document. For all warranty information, including disclaimers, exclusions, terms, conditions and related provisions, refer to the "Arthrex U.S. Product Warranty" section of the Arthrex, Inc. website, found at www.arthrex.com whose provisions are incorporated herein by reference.