English

A. DEVICE DESCRIPTION

The Arthrex Suture Button, BicepsButton™, and Pec Button are comprised of metal buttons which may be provided with suture and pre-loaded onto an inserter. The RetroButton® is comprised of a graft loop and a metal button.

The Distal Biceps Repair Implant System consists of a metal button, suture with a straight needle, button inserter, drill pin and a Tenodesis Screw.

The Pec Repair Implant Delivery Systems consist of metal buttons preloaded with sutures and needles, drill pin and button inserter.

The Implant System for Proximal Tenodesis consists of a metal button, button inserter, suture, drill pin, curved needle with a nitinol loop, and a shoehorn cannula.

B. INDICATIONS

The Suture Button and RetroButton are used for fixation of bone to bone or soft tissue to bone, and are intended as fixation posts, a distribution bridge, or for distributing suture tension over areas of ligament or tendon repair. Specifically, for ACL repair.

- Arthrex recommends the use of #2 o #5 FiberWire® or equivalent with the Suture Button.

- The BicepsButton and Pec Button are used for fixation of bone to bone or soft tissue to bone, and are intended as fixation posts, a distribution bridge, or for distributing suture tension over areas of ligament or tendon repair in the knee, shoulder, and elbow and may include the following indications: anterior cruciate ligament, posterior cruciate ligament, pectoralis repair (minor/major), biceps tendon repair and reattachment (distal/proximal), acromioclavicular repair, and ulnar collateral ligament reconstruction.
- Arthrex recommends the use of FiberTape® with the Large Pec Button, and #2 FiberWire or equivalent with the Pec Repair Button, BicepsButton, and Proximal BicepsButton.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any acute infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Allergies and other reactions to device materials.
- Silicone sensitivity, though very rare, has been reported.

E. WARNINGS

- Do not add additional suture to the RetroButton. The extra suture may impede passage of the device through the femur.
- An internal fixation device must never be reused.
- Do not re-sterilize this device. **Note:** Metal buttons which are provided alone can be re-sterilized if their sterility is compromised.
- All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.

F. MRI SAFETY INFORMATION

I. MR CONDITIONAL

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Suture Buttons manufactured from Titanium or Stainless Steel are MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the Suture Buttons are expected to produce a maximum temperature rise of up to 1.7 °C after 15-minutes of continuous scanning.
- artificial Information**

In non-clinical testing, the image artifact caused by the Suture Buttons can extend up to approximately 17 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

II. MR SAFE

Arthrex Suture manufactured from UHMWPE and polyester with or without silicone elastomer coating, cyanoacrylate, and nylon are MR safe.

G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- RetroButton only:** The loop length must be adequate for the button to exit the cortex.
- RetroButton only:** Measurements of intraosseous length and socket depth are critical for proper implant selection.
- Loop and button should be properly oriented during passing and fixation, per technique guide, in order to function properly.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.

M. INFORMATION

Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

I. MANUAL CLEANING AND DISINFECTION
<ol style="list-style-type: none">Immerse the device in an enzymatic or alkaline cleaning detergent solution. Cleaning solutions can include, but are not limited to: ENZOL® enzymatic, neodisher® Mediclean forte, and Thermosept® alka clean. CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinatedethylenepropylene), ABS (Acrylonitrile Butadiene Styrene), Ultem™, Lexan™, and Cycloc™. If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device. Cleaning solutions should always be mixed to the manufacturer’s specification for concentration and cleaning shall be conducted at ambient temperature unless otherwise stated in the cleaning solution manufacturer’s instructions.
<ol style="list-style-type: none">Scrub device with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the instrument. Immerse the instrument in detergent, agitate and allow it to soak for at least one minute. Rinse the instrument thoroughly with cold distilled water for at least one minute following the cleaning process.
<ol style="list-style-type: none">Immerse devices in disinfection solutions for a minimum of 20 minutes. Suitable disinfection solutions can include, but are not limited to: Cl-DIX®, WAVICIDE®-01, Gigasept®, Kohrosolin®, and equivalent products). Use the supplier’s instructions for preparing the solution. CAUTION: Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinatedethylenepropylene), ABS (Acrylonitrile Butadiene Styrene), Ultem™, Lexan™, and Cycloc™. If non-neutral pH disinfection chemistries are utilized, care should be taken to ensure appropriate rinsing and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device. Disinfection solutions should always be mixed to the manufacturer’s specification for concentration.

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