# English

#### A. DEVICE DESCRIPTION

The ECLIPSE™ Shoulder Prosthesis is a shaft-free humeral joint (hemi-shoulder) prosthesis that is designed as a humeral head replacement device. The shaft-free humeral joint implant has three design elements 1) a humeral head; 2) a trunion; and 3) a hollow screw.

The ECLIPSE Shoulder Prosthesis, when used as a total shoulder joint replacement device, is mated with an Arthrex glenoid implant.

#### B. INDICATIONS

The Arthrex ECLIPSE Shoulder Prosthesis is indicated as a total or hemi shoulder joint replacement where the humeral head and neck are of sufficient bone stock and the rotator cuff is intact or re-constructible for the following specific indications:

- Degenerative pathologies: osteoarthritis, arthritis due to instability, and idiopathic avascular humeral head necrosis
- Post-traumatic pathologies: post-traumatic arthritis, post-traumatic instability, and post-traumatic humeral head necrosis

The Arthrex ECLIPSE Shoulder Prosthesis is intended for uncemented use only.

**Caution:** Federal law restricts this device to sale by or on the order of a physician

#### C. CONTRAINDICATIONS

- Insufficient quantities or quality of humeral head and/or humeral neck bone stock.
- A rotator cuff that is not intact and irreparable.
- Irreducible 3- and 4-part proximal humeral fractures.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period, including severe neuro-arthropathy.
- Do not use for surgeries other than those indicated.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.

#### D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Allergjes and other reactions to device materials.
- Loosening of the implant as a result of changed conditions in load transfer, respectively, fatigue wear or tissue reaction to the implant. Loosening is frequently a consequence of one or several of the above listed risk conditions, but can also be caused by inadequate anchoring technique (see below).
- Dislocation, subluxation, or inadequate scope of movement as a result of failure to achieve optimum positioning of the implant.
- Bone fractures as a result of one-sided overload or weakened bone structure.
- Temporary or permanent nerve damage as a result of pressure or hematoma.
- Cardiovascular complications including venous thrombosis, pulmonary embolism, and cardiac arrest.
- Wound hematoma and delayed wound healing.

#### E. WARNINGS

- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Detailed instructions on the use and limitations of the device, the patient leaflet and the patient implant card should be given to the patient.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- The following operative situations may cause premature loosening and complications:
  - Extreme weakening of the bone structure in preparing the bone bed;
  - Unsuitable selection of the implant size;
  - Inadequate cleaning of the bone bed prior to implantation; and,
  - Excessive use of force in placing or fastening the implant, provoking splintering fractures, or causing the bone to tear.
- An internal fixation device must never be reused.
- Do not re-sterilize this device.

- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- The appropriate Arthrex delivery system is required for proper insertion of the implant.
- Only Arthrex delivery systems, instruments, and trial prostheses should be used for the implantation procedure.
- Endoprostheses may not be processed mechanically or changed in any other way.
- Do not implant any parts that have been scratched or damaged.
- An artificial joint is subject to wear and/or can loosen over a period of time. Wear and loosening may make it necessary to re-operate on an artificial joint.
- An infection in an artificial joint may lead to implant removal.
- This device should only be used in conjunction with other implants designed specifically for use with this system.
- TPS (titanium plasma sprayed)/CaP coated device - Fluid contact other than patient’s blood should be avoided to achieve best ongrowth results.
- Follow your institutions policy for safe disposal of all needles and other sharps or medical waste.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

#### F. MRI SAFETY INFORMATION

##### 1. MR Conditional

*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex Eclipse Shoulder Prosthesis is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the Arthrex Eclipse Shoulder Prosthesis is expected to produce a maximum temperature rise of 3 °C after 15-minutes of continuous scanning.

##### a. Artifact Information

*In non-clinical testing, the image artifact caused by the Arthrex Eclipse Shoulder Prosthesis extends approximately 60 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

#### G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

#### H. FACTORS AND RISKS IMPACTING THE SAFETY AND SERVICE LIFE OF THE IMPLANT

- Patient weight. An overweight patient may present additional risk.
- Extreme stress or strain resulting from work or sport related activity.
- Patients with increased risk of fractures due to repeated strain or trauma, or medical conditions that increase the patient’s risk for trauma, including falls.
- Osteoporosis or osteomalacia.
- Exposure to infectious diseases with possible manifestation in the joints.
- Deformation of the operative site, which can prevent or impede anchoring of the implant.
- Tumors that weak the supportive structure.
- Allergic reactions to implant materials.
- Thromboses and pulmonary infarction caused during implant site preparation and subsequent /treatment.

#### I. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- All prosthesis parts and components should be stored unopened, in the original packaging, and should be kept in protective packaging until immediately prior to use.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolsglossary](http://www.arthrex.com/symbolsglossary).

#### J. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

#### K. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The device consists of a hollow screw, and a trunion manufactured of titanium alloy and a humeral head manufactured of cobalt chromium (Co-Cr) alloy.

Trunions are available either uncoated, or with a TPS/ CaP coating.

**Worldwide except US and Canada:** Humeral Heads are available either uncoated, or with a TiN/N coating.

#### L. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date. Joint prostheses must be treated with great care, prior to implantation. Scratches or nicks on the joint surface of the prosthesis

can result in excessive warp and complications. Do not remove the device from the package until immediately before use.

#### M. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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