English
<div><div><div><div><div><div></div></div></div></div></div></div>
<div><div><div><div><div></div></div></div></div></div>

A. DEVICE DESCRIPTION

The Universal Glenoid™ is a convertible metal-backed glenoid prosthesis which may be used with a UHMWPE inlay bearing to replace the natural humeral glenoid surface, or with a glenosphere for reverse shoulder application. The Universal Glenoid is available in three sizes and is fixeded with central and peripheral fixation screws.

The Unvers Revers™ Porous Coated Baseplate and Inlay Bearing consist of a porous coated titanium glenoid baseplate and a UHMWPE inlay bearing. The baseplate is designed to be used as the glenoid component of the existing Unvers Revers™ Shoulder Prosthesis System or Unvers™ II Shoulder Arthroplasty System.

B. INDICATIONS

The Universal Glenoid is indicated for use in reverse total shoulder arthroplasty in grossly rotator cuff deficient glenohumeral joint with severe arthropathy or a previously failed joint replacement with a gross rotator cuff deficiency. The patient’s joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device.

The Universal Glenoid is indicated for primary, fracture, or revision total shoulder replacement for the relief of pain and significant disability due to gross rotator cuff deficiency, and is intended to be used with the Unvers Revers Shoulder Prosthesis System.

The Universal Glenoid baseplate is Cap (Calcium Phosphate) coated and is intended for cementless use with the addition of screws for fixation.

Outside the US only: The Universal Glenoid is indicated for use in anatomic joint replacement(s) when conditions including severe pain or significant disability resulting from degenerative, rheumatoid, traumatic disease, fracture or injury of the glenohumeral joint; non-union humeral head fractures of long duration; avascular necrosis of the humeral head; neoplastic or dysplastic diseases; or other difficult clinical management problems where arthrodesis or resectional arthroplasty is not acceptable.

For anatomic joint replacement, the Universal Glenoid is indicated to be used as a hemi shoulder replacement, or with the humeral components of the Unvers II or ECLIPSE™ system for total shoulder arthroplasty. The Unvers Revers Porous Coated Baseplate and Inlay Bearing is designed to be used as the glenoid component of the existing Unvers Revers Shoulder Prosthesis System or Unvers II Shoulder Arthroplasty System.

The Unvers Revers Porous Coated Baseplate and Inlay Bearing is indicated for use in a grossly rotator cuff deficient glenohumeral joint with severe arthropathy or a previously failed joint replacement with a gross rotator cuff deficiency. The patient’s joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device.

The Unvers Revers Porous Coated Baseplate and Inlay Bearing is indicated for primary, fracture, or revision total shoulder replacement for the relief of pain and significant disability due to gross rotator cuff deficiency.

The Unvers Revers Porous Coated Baseplate and Inlay Bearing is porous coated and is intended for cementless use with the addition of screws for fixation.

The Unvers Revers Porous Coated Baseplate and Inlay Bearing, and Universal Glenoid Inlay are indicated for use in anatomic joint replacement(s) when conditions including severe pain or significant disability resulting from degenerative, rheumatoid, traumatic disease, fracture or injury of the glenohumeral joint; non-union humeral head fractures of long duration; avascular necrosis of the humeral head; neoplastic or dysplastic diseases; or other difficult clinical management problems where arthrodesis or resectional arthroplasty is not acceptable.

For anatomic joint replacement, the Unvers Revers Porous Coated Baseplate and Inlay Bearing, and Universal Glenoid Inlay are indicated to be used with the humeral components of the Arthrex Unvers II Shoulder Prosthesis System for total shoulder arthroplasty. The Universal Glenoid Inlay is intended for use with the Unvers Revers Porous Coated Baseplate.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any acute infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.

D. ADVERSE EFFECTS

- Infections, both deep and superficial. Foreign body reactions.
- Injury to nerves, blood vessels, muscles, tendons.
- Malunion/nonunion of the bone.
- Necrosis of the humeral head.
- Risk of chronic pain.
- Allergies and other reactions to device materials.
- Loosening of the implant as a result of changed conditions in load transfer, respectively, fatigue wear or breakage of the cement bed and/ or tissue reaction to the implant. Loosening is frequently a consequence of one or several of the above listed risk conditions. In addition a result of failure to achieve optimum positioning of the implant can cause dislocation, subluxation, or inadequate scope of movement.
- Dislocation, subluxation, or inadequate scope of movement as a result of failure to achieve optimum positioning of the implant.
- Bone fractures as a result of one-sided overload or weakened bone substance.
- Cardiovascular diseases including venous thrombosis, pulmonary embolism, and cardiac arrest.
- Wound hematoma and delayed wound healing.
- Tissue reactions caused by allergic reactions to the implanted material, particularly metal, or

caused by accumulations of wear particles or cement particles.

E. WARNINGS

- Failure to achieve the appropriate torque requirements when tightening locking screws may result in the premature loosening of the device.
- Postoperatively, until healing is complete the fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the implant.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management. Removal of the device should be performed using standard surgical practices for device removal.
- Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device. The following operative situations may cause premature loosening and complications:
 - Extreme weakening of the bone structure in preparing the bone bed;
 - Unsuitable selection of the implant size;
 - Inadequate cleaning of the bone bed prior to implantation; and,
 - Excessive use of force in placing or fastening the implant, provoking splintering fractures, or causing the bone to tear.
- An internal fixation device must never be reused. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Do not resterilize this device.
- The operation is to be planned based on the pre-operative X-rays.
- The specific Arthrex implantation instruments are to be used both for preparing the bone bed and for adjusting and inserting the joint prosthesis. The appropriate Arthrex delivery system is required for proper insertion of the implant.
- Only Arthrex delivery systems, instruments, and trial prostheses should be used for the implantation procedure.
- Endoprostheses may not be processed mechanically or changed in any other way.
- Do not implant any parts that have been scratched or damaged. Do not implant any parts that have been scratched or damaged. An artificial joint is subject to wear and/or can loosen over a period of time. Wear and loosening may make it necessary to re-operate on an artificial joint.
- An artificial joint is subject to wear and/or can loosen over a period of time. Wear and loosening may make it necessary to re-operate on an artificial joint.
- An infection in an artificial joint may lead to implant removal.
- This device should only be used in conjunction with other implants designed specifically for use with this system.
- Proper anchoring is of decisive importance for firm, permanent positioning of the prosthesis.
- In the case of joint endoprosthesis intended for cemented anchoring, the surgeon must comply with the instructions and recommendations of the cement manufacturer when it comes to preparation and cementing techniques. Failure to properly align and completely seat the components together can lead to disassociation. Proper technique must be followed to ensure there is no bony or soft tissue interference between modular components. All screws must be adequately tightened to ensure they are not proud to prevent a mechanical interference between modular components. Thoroughly clean and dry tapers, prior to attachment of modular components to avoid crevice corrosion and improper seating. Glenosphere forces are required to verify integrity of the Morse Taper connection between Glenosphere and Baseplate.

F. FACTORS AND RISKS IMPACTING THE SAFETY AND SERVICE LIFE OF THE IMPLANT

- Patient weight. An overweight patient may present additional risk.
- Extreme stress or strain resulting from work or sport related activity.
- Patients with an increased risk of fractures due to repeated strain or trauma, or medical conditions that increase the patient’s risk for trauma, including falls.
- Osteoporosis or osteomalacia.
- Exposure to infectious diseases with possible manifestation in the joints.
- Deformation of the operative site, which can prevent or impede anchoring of the implant.
- Tumors that weaken the supportive structure.
- Allergic reactions to implant materials.
- Thromboses and pulmonary infarction caused during implant site preparation and subsequent treatment.

G. PRECAUTIONS

- Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Patient weight. An overweight patient may present additional risk.
- Extreme stress or strain resulting from work or sport related activity.
- Patients with increased risk of fractures due to repeated strain or trauma, or medical conditions that increase the patient’s risk for trauma, including falls.
- Osteoporosis or osteomalacia.
- Exposure to infectious diseases with possible manifestation in the joints.
- Deformation of the operative site, which can prevent or impede anchoring of the implant.
- Tumors that weaken the supportive structure.
- Allergic reactions to implant materials.

- Thromboses and pulmonary infarction caused during implant site preparation and subsequent /treatment.
- If a revision to the Revers prosthesis is performed and the glenosphere is removed from the baseplate for any reason, a new glenosphere should be implanted.
- If a revision to the Revers prosthesis is performed and the baseplate is not removed, the baseplate Morse Taper socket should be inspected for damage. If damage is identified or suspected, the baseplate should be removed and replaced.

H. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex Universal Glenoid is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the Ar-threx Universal Glenoid is expected to produce a maximum temperature rise of 3 ° C after 15-minutes of continuous scanning.

a. Artifact Information

In non-clinical testing, the image artifact caused by the Arthrex Universal Glenoid extends approximately 60 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

1. PACKAGING AND LABELING

- Arthrex implants should be accepted only if the factory packaging and labeling arrive intact.
- All prosthesis parts and components should be stored unopened in the original packaging, and should be kept in protective packaging until immediately prior to use.
- Contact customer service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at **www.arthrex.com/symbols glossary**.

J. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DDU-0023-XX and ANSI/AAMI ST79, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities”.

K. MATERIAL SPECIFICATIONS

Refer to the package label for the materials. The device may be made of:

- Calcium Phosphate (CaP)
- Titanium Alloy (Ti6Al4V)
- Titanium Plasma Spray (TPS)
- Ultra-High Molecular Weight Polyethylene (UHMWPE)

L. STORAGE CONDITIONS

This device must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

M. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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