English

A. DEVICE DESCRIPTION
The Interference Screw family includes the Sofi Screw™, Interference Screw, Tenodesis Screw, RetroScrew™, Round Delta Interference Screw, Delta Tapered Interference Screw, FastThread™ Interference Screw and Bio Cortical™ Screw.

These interference screw products may be threaded or partially threaded, cannulated or non-cannulated.

B. INDICATIONS

These devices are intended to be used for fixation of tissue, including ligament or tendon to bone, or a bone to tendon to bone. See below for specific indications. (Tenodesis Screws are not intended to be used for bone to bone fixation.)

See below for specific indications.

Shoulder: Rotator Cuff Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair (Except in EU)

Foot/Ankle (except Delta Tapered Screw): Lateral Stabilization, Medial Stabilization, Flexor Digitorum Longus, Metatarsal Ligament Repair (Except in EU), Heor Hallucis Longus for Achilles Tendon reconstruction (Except in EU), tendon transfers in the foot and ankle

Knee: Anterior Cruciate Ligament Repair (**except 2.5 mm Tenodesis Screw**), Medial Collateral Ligament Repair (Except in EU), Lateral Collateral Ligament Repair (Except in EU), Iliotibial Band Tenodesis (Except in EU), Posterior Cruciate Ligament Repair (Except in EU), **BioComposite Interference Screw, FastThread Interference Screw, Round Delta BioComposite Interference Screw, Delta Tapered BioComposite Interference Screw, BioComposite RetroScrew**

Elbow (except Delta Tapered screw): Biceps Tendon Reattachment, Ulnar/Radial Collateral Ligament Reconstruction (Except in EU)

Hand/Wrist (except Delta Tapered screw): Carpopecta-rcarpal joint arthroplasty (scapho thumb arthroplasty) (**except 2.5 mm Tenodesis Screw**), Basal thumb Ligament Reconstruction (Except in EU), Ulnar/Radial Collateral Ligament Reconstruction (Except in EU)

C. CONTRAINDICATIONS

Bio-Tenodesis **only:** Screw sizes that are smaller than 7mm may not be appropriate for the laser indication.
Insufficient quality or quality of bone.
Blood supply limitations and previous infections which may retard healing.
Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
Any active infection or blood supply limitations.
Conditions that limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not disturb, distort or disrupt the growth plate.
Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS
1. Infections, both deep and superficial.
2. Foreign body reactions.
3. Allergic-like reactions to PLA materials (PLLA, PLLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.

E. WARNINGS
1. Caution: Federal law restricts this device to sale by or on the order of a physician.
This device is intended to be used by a trained medical professional.
3. An internal fixation device must never be reused.
4. **Metal implants only:** All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
Bioresorbable implants: Do not re-sterilize this device.
Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
8. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.

9. Metal implants only: All metallic implant devices used for this surgical procedure should have the same metallurgical composition.

10. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.

11. Metal Implants Only: Removal of supplemnetal fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to the stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.

12. Follow your institutions policy for safe disposal of all needles and other sharps or medical waste.

13. Biobazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
14. Biohazardous materials should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional
Non-clinical testing and in-vivo electromagnetic simulations indicate that the metal (titanium and stainless steel) Interference Screws are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only.
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less.
- Maximum MR system reported, but within averaged specific absorption rate (SAR) of 1-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system.

Under the scan conditions defined, the Interference Screws are expected to produce a maximum temperature rise of up to 6 °C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Interference Screw can extend up to approximately 120 mm from this implant when imaged using a gradient echo pulse sequence and a 3.0 Tesla MR system.

2. MR Safe
The Interference Screws devices manufactured from only poly-etheretherketone (PEEK), polypropylene or poly (L-Lactide), poly (L-Lactide, PLLDLA), biphasic calcium phosphate, biphasic tricalcium phosphate (TCP) and/or hydroxyapatite (HA) are MR safe.

G. PRECAUTIONS

1. Seek the professional judgement when determining the appropriate screw size based on the specific indication, preferred surgical technique, and patient history.

2. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides

detailed surgical techniques in print, video, and electronic formats. The Arthrex website provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an on-site demonstration.

Under insertion of the device may leave the proximal end of the implant protruding beyond the cortical bone, which could potentially cause soft-tissue irritation or pain post-operatively.
4. Bio-Tenodesis Screw only: Use the appropriate size Arthrex Drill to create a pilot hole in the bone.
5. Bio Cortical and Delta Tapered Interference Screw only: Insert the screwdriver into the screw to a fully seated position. Failure to engage the screw completely may result in damage to the implant and/or cancellous or non-cannulated.
6. Bio-absorbable Interference Screw only: It is important to completely seat the screwdriver to prevent potential stripping of the hex and/or screw fracture during insertion or removal.
7. Bio-absorbable Interference Screw only: If inserting the Interference Screw through the anterior tibial portal, use flexion angle of 120° must be maintained throughout the entire insertion process. **Not maintaining or changing the knee flexion angle during screw insertion may result in screw divergence or failure of the screwdriver.** If achieving and maintaining the appropriate flexion angle is not possible or reasonable, a central transpatellar tendon portal should be considered for proper insertion.
8. Bio-absorbable Interference Screw only: During screw size changes or revision procedures, the sheath may be reinserted with the screwdriver over a guide pin to back out the screw in the case of a jammed screw.

H. PACKAGING AND LABELING
1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolslogos.

I. VALIDATION

The recommended cleaning, disinfection, and sterilization methods in this DFU have been validated in compliance with federal and international guidance/standards. In accordance with ISO 17665, the “overall” approach was used for sterilization validation, and demonstrates a sterility assurance level (SAL) of 10⁻⁶. Cleaning, disinfection, and sterilization equipment and methods vary in performance characteristics. Therefore, it is the responsibility of the facility user to perform the appropriate validation testing for any use beyond recommended performance characteristics.

In accordance with EN ISO 17664 and ANSI/ASTM, limit loads and a means for assessing chemical residuals following cleaning have been established for the device. Arthrex devices should be tested for residuals following the manual cleaning and disinfection process or the machine (automated) cleaning and disinfection process, a clinically relevant method was utilized for testing the safety of residuals as part of the validation protocol. Deionized (critical) water was utilized as the terminal rinse water to ensure that residuals will not interfere with subsequent processing steps.

Repeated processing has minimal effect on these devices. End of life is normally determined by wear and damage due to the intended use. The user assumes liability and is responsible for the use of a damaged and dry device.

A device labeled as a Single Use device must never be reused. Used refers to those single use devices that have come into contact with blood, urine, or other body fluids. Any unused single use device that has been exposed to blood, urine, tissue, or body fluids must not be reprocessed and must be discarded.

The instructions in this DFU were developed using the guidance given in the following standards:

- ANSI/ASTM S779, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities”
- ISO 17664: Sterilization of medical devices – Information for medical device manufacturers for the processing of resterilizable medical devices
- ISO 17665-1: Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices
- ANSI TR80-2011: A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices
- ANSI S777: Containment devices for reusable medical device sterilization

J. CLEANING AND DISINFECTION

Certain Arthrex devices that may be used during a procedure are cleaned non-sterile and must be adequately cleaned and sterilized prior to use or re-use. All devices are to be cleaned, disinfected, and sterilized prior to each application; this is required as well for the first use after delivery of the unsterile devices. An effective cleaning and disinfection is an indispensable requirement for an effective sterilization of the device. Single-use devices must be cleaned separately from solid devices.

If possible, the Machine (Automated) procedure should be used for cleaning and disinfection of the instruments. The Manual Cleaning procedure should only be used if an automated procedure is not available; in this case, the significantly lower efficiency and reproducibility of the manual procedure should be considered. The Preliminary Cleaning steps are to be performed in both cases. Manual cleaning may require on-site validation by the healthcare facility and appropriate procedures should be in place to avoid human factor variability.

I. DETERGENT SELECTION
Consider the following points during selection of the cleaning detergent:

1. Suitability of the cleaning agent for ultrasonic cleaning (no foam development).
2. Compatibility of the cleaning agent with the instruments. Arthrex recommends the use of neutral pH or enzymatic cleaning agents. Alkaline agents may be used to clean devices in countries where required by law or local ordinance, or where prior diseases such as Transmissible Spongiform Encephalopathy (TSE) or Creutzfeldt – Jakob disease (CJD) are a concern. **Caution: Low acid or alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluoroethylene Propylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Lexan™. If use neutral pH cleaning chemistry are utilized, care should be taken to ensure appropriate rinsing, as validated by the end-user facility, and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.**
Follow the instructions of the detergent manufacturer regarding use concentration and temperature for either manual or automated cleaning. Use only freshly prepared solutions as well as only purified/ highly purified water as both for final rinse, and a soft, low-linting cloth or filtered medical grade air for drying, respectively.

II. PRELIMINARY CLEANING
Note: Do not assemble/disassemble these devices is required unless stated on the labeling, directions for use, or literature assembly instructions (LAI) pertaining to cleaning, disinfection, and sterilization. Devices that require disassembly are to be disassembled prior to cleaning.

1. Remove excess soil from devices, especially in areas such as joints and crevices, by cleaning the surfaces with a sponge or brush under cold running water or with a non-shedding disposable wipe for a minimum of 30 seconds.

2. Rinse the devices at least 1 minute under running utility water (temperature <= 35 °C/95 °F). Special attention should be given to handles, joints, crevices, and other hard-to-reach areas.
3. Immerse the devices in cleaning solution using an ultrasonic bath. While immersed in solution, brush the devices for 1 minute using a soft-bristled brush. Special attention for 1 minute using a soft-bristled brush. Special attention for 1 minute using a soft-bristled brush. Special attention for 1 minute using a soft-bristled brush.

4. After brushing, turn on ultrasonic power and soak and sonicate for 10 minutes at a minimum of 40± 5 Hz. Ensure contact with the device in the open position and that lumens have complete contact with cleaning solution during soaking.
5. Remove the devices from the cleaning solution and rinse at least 1 minute with utility water. Thoroughly and aggressively rinse lumens, joints, crevices, and other hard-to-reach areas. Complete the cleaning, and drying of the end user has the option to perform either Manual Cleaning and Disinfection or Machine (Automated) Cleaning and Thermal Disinfection (preferred).

III. MACHINE (AUTOMATED) CLEANING AND DRYING

Considerations for the selection of the washer-disinfector:

- Capable of providing an approved program for thermal disinfection (appropriate exposure time and temperature according to A concept)
- Final rinse completed with purified (critical, eg. RO or DI) water, and utilizes only filtered air for drying.

After preliminary cleaning is complete, load the devices in the washer-disinfector such that all design features of the device are accessible to cleaning and such that design features that might retain liquid can drain (for example, hinges should be open and cannulations/holes positioned to drain).
2. If using alkaline cleaning agents, a neutralization step should be utilized as appropriate.
3. Follow the manufacturer's cycle with fundamentally approved efficiency of the washer-disinfector (for example, CE marking according to EN ISO 15885 or FDA approval/clearance/registration). The following minimum recommended automatic wash cycle parameters were utilized by Arthrex during the validation of these instruments:

RECOMMENDED WASHING CYCLE PARAMETERS				
Phase	Recirculation Time	Temperature	Detergent	
Pre-Wash	3 Minutes	Cold Water	N/A	
Cleaning Wash	10 Minutes	Follow detergent manufacturer's recommendation	Enzymatic or alkaline detergent	
Neutralization Rinse (optional)	2 Minutes	Follow detergent manufacturer's recommendation	Neutralizing agent (as needed)	
Rinse	3 Minutes	Cold Water	N/A	
Thermal Disinfection Rinse	5 Minutes	90°C (194°F)	N/A	
Drying	Minimum 6 Minutes or until visibly dry	Minimum 100°C (212°F)	N/A	

4. Remove the devices from the washer-disinfector following the completion of the program and check devices for visible soil. Reinspect if any soil is visible and re-inspect; otherwise, proceed to sterilization.

IV. MANUAL CLEANING AND DISINFECTION

Following preliminary cleaning, the instructions for Manual Cleaning and Disinfection may be followed as an alternative cleaning method to Machine (Automated) Cleaning and Thermal Disinfection if an automated procedure is not available.

1. After preliminary cleaning is complete, repeat steps 1-5 in the Preliminary Cleaning section within this DFU, including rinsing, immersion and sonication, and post-rinsing. Final rinsing should be completed with purified (critical, eg. RO or DI) water.
2. Check devices for visible soil. Repeat cleaning if any soil is visible and re-inspect.
3. Soak the devices for the given soaking time (provided by the disinfectant manufacturer) in disinfectant solution so that the devices are sufficiently covered. Make sure that there is no contact between the devices. Ensure that the device is in the open position during soaking. Attach movable parts at least five times during disinfection, as applicable.

4. Remove the devices from the disinfectant solution and rinse per disinfectant manufacturer's instructions.

5. Dry devices thoroughly using filtered medical grade air or a soft, low-linting cloth. Proceed to Sterilization section.

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M. STORAGE CONDITIONS
Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.
Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of sterile devices is not limited; the devices are manufactured from non-sterile materials which do not raise any question of device stability when stored under recommended conditions.

It is the responsibility of the end-user to ensure devices, once sterilized, are stored in such a way as to maintain the sterility of the device until use. Sterile, packaged devices should be stored in the limited access area that is well ventilated and provides protection from

dust, moisture, insects, and temperature/humidity extremes. Sterile devices packages should be carefully examined prior to opening and that package integrity has not been compromised. Maintenance of sterile package integrity is generally event related. If a sterile wrap is torn, show any evidence of tampering, or has been exposed to moisture, the device or part must be cleaned, repackaged, and sterilized.

N. INFORMATION

1. **In EU Only:** Procedures carried out using these devices may be subject to the specific national regulatory requirements for such devices.
In EU Only: The clinical benefits associated with the use of these devices outweigh the known clinical risks.
3. **In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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