# English

**A. DEVICE DESCRIPTION**  
The Spiked and Spikeless Ligament Staps are metal devices. The Small Staps are metal devices.

The DynaNite® Nitinol Staple and DynaNite® SuperMX Staple provide a means of bone fixation for the management of fractures and reconstructive surgery.

- DynaNite® Staps are designed to exhibit super-elastic properties.
- The Staps are sold in two configurations:
  - as part of a single use disposable kit – the kit includes a staple, pre-loaded onto a delivery device, a drill bit, a drill guide, alignment pins and a tamping and
  - as a staple loaded on a delivery device – for procedures requiring additional same sized staples.

**B. INDICATIONS**  
The DynaNite Nitinol Staple, DynaNite SuperMX Staple and Small Staps are intended to be used for fixation such as: Lisfranc arthrodesis, mono- or bi-cortical osteotomies in the forefoot, first metatarsals, akin osteotomy, midfoot and hindfoot arthrodeses or osteotomies, fixation of osteotomies for hallux valgus treatment (Scarf and Chevron), and arthrodesis of the metatarsometatarsal joint to reposition and stabilize metatarsus primus varus.

The Spiked and Spikeless Ligament Staps are used in the repair of ligamentous tissue (only the ALG).

**C. CONTRAINDICATIONS**

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.
- DynaNite Staple only:** Contaminated bone surface that would not allow for staple placement.
- DynaNite Staple only:** Pathologic bone conditions such as osteopenia that would impair the ability to securely fix the implant.

**D. ADVERSE EFFECTS**

- Infections, both deep and superficial.
- Foreign body reactions.

**E. WARNINGS**

- Caution: Federal law restricts this device to sale by or on the order of a physician
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be reused.
- All metal implant devices used for this surgical procedure should have the same metallurgical composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system and procedure for proper implantation of the device should be followed. Any decision to know the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- If any of the following conditions are present—nomion, osteoporosis, a markedly unstable comminuted fracture or any of the factors listed in the Contraindications and/or Warnings—Precautions sections, then the following can occur: loosening, bending, cracking, fracture of the staple or loss of fixation in bone.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management and avoid re-fracture.
- Detailed instructions on the use and limitations of this device, the patient leaflet and the patient implant card should be provided to the patient.
- DynaNite Staple only:** The implants are not designed to replace normal healthy bone or withstand the stress placed upon the device by full or partial weight-bearing or load bearing in the presence of nonunion, delayed union or incomplete healing. Immobilization of the treatment site using routine methods (casting, splinting, etc.) should be maintained until bone healing has occurred (4-6 weeks).
- DynaNite Staple only:** Reduction of the fracture should be achieved and maintained prior to activating the device. The staple should not be relied upon to achieve closure or reduction of the fracture line.
- DynaNite Staple only:** Additional processing or reprocessing of the implant, may affect the material properties of the Nitinol, changing or otherwise reducing the effectiveness of the implant.
- DynaNite Staple only:** Drill bits are designed for use specifically with the kits supplied drill guides. Reprocessing may affect the compatibility with other instruments and usability of the instruments.
- DynaNite Staple only:** If sterilization is compromised prior to insertion, a different sterile implant will not be used. Do not re-sterilize or re-use explained devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

**F. MRI SAFETY INFORMATION**

*Non-clinical testing and in-vitro electro-magnetic simulations demonstrated that the Arthrex Staps are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
  - Magnetic spatial gradient magnetic field of 3000 Gauss/cm or less
  - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 1-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the Arthrex Staps is expected to provide a maximum temperature rise of up to 6°C after 15-minutes of continuous scanning.
- In non-clinical testing, the image artifact caused by the Arthrex Staps can extend up to approximately 120 mm from the implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*
- G. PRECAUTIONS**

Surgons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or contact your Arthrex representative for an onsite demonstration.

2. Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.

- For proper fixation, an adequate amount of bone must be between the device and the fracture site.
- For hard cortical bone, select the appropriate size guidewire or drill bit for the device.

**H. PACKAGING AND LABELING**

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolslogosary](http://www.arthrex.com/symbolslogosary).

**I. VALIDATION**

The recommended cleaning, disinfection and sterilization methods in this DFU have been validated in compliance with federal and international guidance standards. In accordance with ISO 17665, the "sterilize" approach to validation, and validation, and demonstrates a sterility assurance level (SAL) of 10<sup>-6</sup>. Cleaning, disinfecting, and sterilizing equipment and materials vary in performance characteristics. Therefore, it is the responsibility of the facility/ end user to perform the appropriate validation testing for any use beyond recommended performance characteristics.

In accordance with EN ISO 17664 and ANSI TR190, limits values and a means for monitoring chemical residues following cleaning have been established for the product. In assessing the level of cleaning residuals following the manual cleaning and disinfection process or the machine (automated) cleaning and disinfection process, a clinically relevant method was utilized for testing the safety of residuals as part of the validation protocol. Distilled (critical) water was utilized as the terminal rinse water quality to ensure that residuals will not interfere with subsequent processing steps.

Repeating processing has minimal effect on these devices. End of use is normally determined by wear and damage due to the intended use. The user assumes full liability and is responsible for the use of a damaged and/or deformed device.

A device labeled as a Single Use device must never be reused. Users refers to those single use devices that have come into contact with blood, bone, tissue, or other body fluids. Any unused single use device that has been exposed to blood, bone, tissue, or body fluids must not be re-sterilized or re-used.

The instructions in this DFU were developed using the guidance given in the following standards:

- ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities"
- ISO 17664, Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices
- ISO 17665-1: Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices
- AMI TR18-2011: A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices
- AMI ST77: Containment devices for reusable medical device sterilization

**J. CLEANING AND DISINFECTION**

Certain Arthrex instruments that are used during this procedure are cleaned non-sterile and must be adequately cleaned and sterilized prior to re-use or re-use. All devices are to be cleaned, disinfected, and sterilized prior to each application; this is required as well for the first use after delivery of the unsterile devices. An effective cleaning and disinfection is an indispensable requirement for an effective reprocessing of the device. Single-use devices must be cleaned separately from solid devices.

If possible, the Machine (Automated) procedure should be used for cleaning and disinfection of the instruments. The Manual Cleaning procedure should only be used if an automated procedure is not available; in this case, the significantly lower efficiency and reproducibility of the manual procedure should be considered. The Preliminary Cleaning steps are to be performed in both cases. Manual cleaning may require onsite validation by the healthcare facility and appropriate procedures should be in place to avoid human factor variability.

**DETERGENT SELECTION**  
Consider the following points during selection of the cleaning detergent:

- Suitability of the cleaning agent for ultrasonic cleaning (no foam development).
- Compatibility of the cleaning agent with the instruments. Arthrex recommends the use of neutral pH or enzymatic cleaning agents. Alkaline agents may be used to clean devices in countries where required by law or local ordinance, or where prior diseases such as Transmissible Spongiform Encephalopathy (TSE) or Creutzfeldt - Jakob disease (CJD) are a concern. **Caution: Low acid or alkaline solutions are not recommended, as they corrode metal parts and modified aluminum and compromise polymer plastics, such as PEP (Fluorinatedethylene/propylene), ABS (Acrylonitrile Butadiene Styrene), Ultem™, Lexan™, and Cyocel™. If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing as validated by the end-user facility, and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.**
- Follow the instructions of the detergent manufacturer regarding use concentration and temperature for either manual or automated cleaning. Use only freshly prepared solutions as well as only purified, highly purified water at least for final rinse, and a soft, low-linting cloth or/or filtered medical grade air for drying, respectively.

**II. PRELIMINARY CLEANING**  
**Note:** No assembly/disassembly of these devices is required unless stated on the labeling, directions for use, or literature assembly in the presence of nonunion, delayed union or incomplete healing. Immobilization of the treatment site using routine methods (casting, splinting, etc.) should be maintained until bone healing has occurred (4-6 weeks).

- Remove excess soil from devices, especially in areas such as joints and crevices, by cleaning the surfaces with a sponge or brush under cold running water or with a non-shedding disposable wipe for a minimum of 30 seconds.
- Rinse the devices at least 1 minute under running utility water (temperature < 35 °C/95 °F). Special attention should be given to lumens, joints, crevices, and other hard-to-reach areas.
- Immense the devices in cleaning solution inside an ultrasonic bath. While immersed in solution, brush the device for 1 minute using a soft-bristled brush. Special attention should be given to lumens, joints, crevices, and other hard-to-reach areas. Lengths should be brushed with appropriate diameter and long bristle sizes for the particular lumen. Actuate movable parts at least (5) times during soaking, as applicable.
- After brushing, turn on ultrasonic power and soak and sonicate for 10 minutes at a minimum of 40 ± 5 kHz. Ensure devices are in the open solution and that lumens have complete contact with cleaning solution during soaking.
- Remove the devices from the cleaning solution and rinse at least 1 minute with utility water. Thoroughly and aggressively rinse lumens, joints, crevices, and other hard-to-reach areas.
- After the completion of preliminary cleaning, the end user has the option to perform either Manual Cleaning and Disinfection or Machine (Automated) Cleaning and Thermal Disinfection (preferred).

**III. MACHINE (AUTOMATED) CLEANING AND THERMAL DISINFECTION**

Considerations for the selection of the washer-disinfector:

- Capable of providing an approved program for thermal disinfection (appropriate exposure time and temperature according to A. compendium).
- Final rinse completed with purified (critical, e.g. RO or DI) water, and utilizes only filtered air for drying
- After preliminary cleaning is complete, load the devices in the washer-disinfector such that all design features of the device are accessible to cleaning and such that design features that might retain liquid can drain (for example, hinges should be open and cannulations/holes positioned to drain).
- If using alkaline cleaning agents, a neutralization step should be utilized as appropriate.
- Run an automated wash cycle with fundamentally approved cleaning agents and disinfectant (for example, ClE marking according to EN ISO 15858 or FDA approved clearance/registration). The following minimum recommended automated wash cycle parameters were utilized by Arthrex during the validation of these instruments.

cycle parameters were utilized by Arthrex during the validation of these instruments.

RECOMMENDED WASHING CYCLE PARAMETERS			
Phase	Recirculation Time	Temperature	Detergent
Pre-Wash	3 Minutes	Cold Water	N/A
Cleaning Wash	10 Minutes	Follow detergent manufacturer’s recommendation	Enzymatic or alkaline detergent
Neutralization Rinse (optional)	2 Minutes	Follow detergent manufacturer’s recommendation	Neutralizing agent (as needed)
Rinse	3 Minutes	Cold Water	N/A
Thermal Disinfection Rinse	5 Minutes	90°C (194°F)	N/A
Drying	Minimum 6 Minutes or until visibly dry	Minimum 100°C (212°F)	N/A

- Remove the devices from the washer-disinfector following the completion of the program and check devices for visible soil. Repeat cleaning if soil is visible and re-inspect; otherwise, proceed to Sterilization section.

**IV. MANUAL CLEANING AND DISINFECTION**

Following Preliminary Cleaning, the instructions for Manual Cleaning and Disinfection may be followed as an alternative cleaning method to Machine (Automated) Cleaning and Thermal Disinfection if an automated procedure is not available.

- After preliminary cleaning is complete, repeat steps 1-5, including in the Preliminary Cleaning section with the following procedure:
  - rinsing, immersion and sonication, and post-rinsing. Final rinsing should be completed with purified (critical, e.g. RO or DI) water.
  - Check devices for visible soil. Repeat cleaning if soil is visible and re-inspect.
  - Soak the devices for the given soaking time (provided by the disinfectant manufacturer) in disinfectant solution so that the devices are sufficiently covered. Make sure that there is no contact between the devices. Ensure that the device is in the open position during soaking. Activate movable parts at least five times during disinfection, as applicable.
  - Remove the devices from the disinfectant solution and rinse per disinfectant manufacturer’s instructions.
  - Dry devices thoroughly utilizing filtered medical grade air or a soft, low-linting cloth. Proceed to Sterilization section.

**K. STERILIZATION**  
This device may be provided either sterile or non-sterile. Check the package labeling for more information. For devices that are not provided in a terminally sterilized configuration, sterilization is to be performed following cleaning, disinfection, and sterile packaging prior to use, and may be re-sterilized (if unused) following cleaning, disinfection, and sterile packaging prior to use.

Certain that are provided in a terminally sterilized configuration should never be re-sterilized under any conditions. Repeat Arthrex instructions that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to re-use. Please refer to DFI-0023-XX and ANSI/AAMI ST79 for specific information.

**I. STERILE PACKAGING**  
Single-Use devices should be packed so as to ensure that the pack is large enough to contain the device without stressing the seals. Packaging should be completed utilizing a pouch or wrap which conforms to the recommended specifications for steam sterilization as outlined below. If a wrap is utilized, it should be completed following AAMI double-wrap or equivalent guidelines with an appropriate wrap. An appropriate wrap is one that, for example, is cleared by the FDA or the local governing body at the point of use. Devices may also be placed into an approved reusable rigid sterilization container. Aesculap SteriContainer™ rigid containers with perforated bottoms and lids are approved for use with Arthrex, Inc. devices.

Sets: Where appropriate, cleaned, disinfected and inspected devices should be placed into trays/cases as provided or in general-purpose sterilization trays. The total weight of trays/cases should not exceed 11.4kg/25 lbs. (other local limits below 11.4kg/25 lbs. may apply). Trays/cases should be double wrapped following AAMI or equivalent guidelines with an appropriate wrap. An appropriate wrap is one that, for example, is cleared by the FDA or the local governing body at the point of use. Sets may also be placed into an approved reusable rigid sterilization container. Aesculap SteriContainer™ rigid containers with perforated bottoms and lids are approved for use with Arthrex, Inc. devices.

Areas designated for specific devices shall contain only devices intended for those areas. Devices should not be stacked or placed in close contact. Only Arthrex devices should be included in the trays or cases. These validated reprocessing instructions are not applicable to trays or cases that include devices not intended to be used with Arthrex trays or cases.

**II. STERILIZATION**  
Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in the table below. Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer’s instructions.

RECOMMENDED STEAM STERILIZATION PARAMETERS					
Cycle Type	Minimum Exposure Temperature	Minimum Exposure Time	Minimum Drying Time <sup>1</sup>	Minimum Cooling Time <sup>2</sup>	
US Pre-vacuum Cycle	132°C (270°F)	4 Minutes	30 Minutes	30 Minutes	
UK Pre-vacuum Cycle	134°C (273°F)	3 Minutes	30 Minutes	30 Minutes	
Prevacuum Cycle (Priority)	134°C (273°F)	18 Minutes	30 Minutes	30 Minutes	

<sup>1</sup>Drying times vary according to load size and should be increased for larger loads.

<sup>2</sup>Cooling times vary according to the sterilizer used, device design, temperature and humidity of ambient environment, and type of packaging used. Cooling process should comply with ANSI/AAMI ST79.

<sup>3</sup>Reprocessing parameters recommended by the World Health Organization (WHO) where there is concern regarding TSE/CJD contamination.

**L. MATERIAL SPECIFICATIONS**  
Refer to the package label for the materials. This device is made of stainless steel, titanium, Cobalt Chromium or Nitinol.

**M. STORAGE CONDITIONS**  
Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise a question of device stability when stored under recommended conditions.

It is the responsibility of the end-user to ensure devices, once sterilized, are stored in such a way as to maintain the sterility of the device until use. Sterile, packaged devices should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, and temperature/humidity extremes. Sterile device packages should be carefully examined prior to opening to ensure that package integrity has not been compromised. Maintenance of sterile package integrity is generally not expected. If a sterile wrap is torn, shows any evidence of tampering, or has been exposed to moisture, the device or set must be cleaned, repackaged, and sterilized.

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