A. DEVICE DESCRIPTION

The BioProWick Postoperative Antimicrobial Wound Dressing and Cold Therapy System contains one cover dressing, a set of five foam pads, two compression straps, one cold pack, and one set of five JumpStart® antimicrobial wound dressings for use during the first 48 hours after surgery. The single size BioProWick Dressing is designed to accommodate most knees and shoulders, and is intended for absorbing body fluids that may exude from incisions after surgery.

B. INDICATIONS

The BioProWick Antimicrobial Wound Dressing is indicated for the management of post-surgical wounds. The BioProWick Cold Therapy Wrap is intended for providing cold therapy during the first 48 hours after arthroscopic or mini-open knee or shoulder surgery.

See the “JumpStart® Contact Layer Antimicrobial Wound Dressing Instructions For Prescription Use”. The Directions/Instructions For Use should accompany the respective products about the location they should be applied.

C. CONTRAINDICATIONS

1. Do not use for surgeries other than those indicated.
2. See the “JumpStart Contact Layer Antimicrobial Wound Dressing Instructions For Prescription Use”. The Directions/Instructions For Use should accompany the respective products about the location they should be applied.

D. ADVERSE EFFECTS

1. See the “JumpStart Contact Layer Antimicrobial Wound Dressing Instructions For Prescription Use”. The Directions/Instructions For Use should accompany the respective products about the location they should be applied.

E. WARNINGS

1. Do not re-sterilize this device.
2. Detailed instructions on the use and limitations of this device should be given to the patient.
3. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.

F. MRI SAFETY INFORMATION

The BioProWick Postoperative Antimicrobial Wound Dressing and Cold Therapy System is MR Unsafe.

G. PRECAUTIONS

1. The BioProWick Dressing should be removed within 48 hours after placement.
2. Do not immerse the dressing in fluid or expose it to excess moisture, as the dressing could become saturated.
3. BioProWick Postoperative Antimicrobial Wound Dressing System is for SINGLE USE ONLY. DO NOT RESTERILIZE OR REUSE.
4. The BioProWick Cold Therapy Wrap is not intended to provide heat therapy. Do not heat the Cold Pack in a microwave oven or by any other method.
5. The BioProWick Cold Therapy Wrap is to be used after application of the BioProWick Postoperative Antimicrobial Wound Dressing.
6. Do not allow the BioProWick Cold Therapy Wrap to get wet.
7. Do not place the Cold Pack directly against skin.

8. The cold therapy should be applied per the hospital/physician protocol.
9. See the “BioProWick Cold Therapy Wrap Directions For Use” and “JumpStart Contact Layer Antimicrobial Wound Dressing Instructions For Prescription Use”. The Directions/Instructions For Use should accompany the respective products about the location they should be applied.
10. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

H. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials. This device is made of:

- The foam pad device consists of Polyethylene, Acrylic Adhesive and Polyurethane.
- The dressing device consists of Polyethylene, Polypropylene, Acrylic Adhesive, Nylon and SuperCore®.
- The cold pack device consists of Polypropylene, Ethyl Vinyl Acetate, Acrylic Adhesive and Nylon.
- The compression strap consists of Polypropylene, Acrylic Adhesive and Nylon.

These products are latex free.

See the “BioProWick Cold Therapy Wrap Directions For Use” and the “JumpStart Contact Layer Antimicrobial Wound Dressing Instructions For Prescription Use”. The Directions/Instructions For Use should accompany the respective products about the location they should be applied.

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

Store in a cool, dry place and in a manner that protects the integrity of the package and sterile barrier.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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