<p align="center">Symbols Used on Labeling Auf der Verpackung verwendete Symbole Symboles utilisés sur l'étiquetage Simboli usati sull'etichetta Símbolos utilizados en las etiquetas Símbolos utilizados na embalagem</p>
--

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A. DEVICE DESCRIPTION

The Arthrex Mixing and Delivery System consists of a piston syringe with a movable plunger and cap to facilitate mixing and delivery; and may include straight and curved (tuohy) delivery needles; a mating obturator for delivery needles; luer connectors; and a funnel to facilitate filling of the syringe barrel. The system will be offered with either a 3 mL or 14 mL syringe barrel and may be provided either empty or pre-filled with allograft, autograft, or synthetic bone graft materials.

B. INDICATIONS

The Arthrex Mixing and Delivery System is indicated for the delivery of allograft, autograft, or synthetic bone graft materials to all orthopedic surgical sites. In addition, it is designed to facilitate pre-mixing of allograft, autograft, or synthetic bone graft materials with I.V. fluids, blood, plasma, platelet rich plasma, bone marrow or other specific blood component(s) as deemed necessary by the clinical use requirements.

C. CONTRAINDICATIONS

1. Blood supply limitations and previous infections, which may retard healing.
2. Any active infection or blood supply limitations.

D. ADVERSE EFFECTS

1. Allergies and other reactions to device materials.

E. WARNINGS

1. Do not re-sterilize this device.
2. Do not reuse this device.
3. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

F. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

G. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.

H. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

I. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

- Needle: Surgical-grade stainless steel. ABS Luer Lock
- Syringe: Polyamide
- Cap: Polypropylene
- Luer Connectors: ABS
- Obturator: Surgical-grade stainless steel, ABS Luer Lock
- Funnel: Polypropylene

J. STORAGE CONDITIONS

This device is supplied sterile. Check the package labeling for the sterilization method. Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

K. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

Deutsch	DFU-0187r3
---------	------------

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