

# iBalance® TKA

## Implant Sizing Compatibility Matrix



# TIBIA

# FEMUR

SIZE		1	2	3	4	5	6	7	8	9	10	
	PART NO.	AR-513-T1	AR-513-T2	AR-513-T3	AR-513-T4	AR-513-T5	AR-513-T6	AR-513-T7	AR-513-T8	AR-513-T9	AR-513-T10	
PS Poly w/PS Femur	1	AR-516-1	PS-1	PS-2	PS-3							
	2	AR-516-2	PS-1	PS-2	PS-3	PS-4						
	3	AR-516-3	PS-1	PS-2	PS-3	PS-4	PS-5					
	4	AR-516-4		PS-2	PS-3	PS-4	PS-5	PS-6				
	5	AR-516-5			PS-3	PS-4	PS-5	PS-6	PS-7			
	6	AR-516-6				PS-4	PS-5	PS-6	PS-7	PS-8		
	7	AR-516-7					PS-5	PS-6	PS-7	PS-8	PS-9	
	8	AR-516-8						PS-6	PS-7	PS-8	PS-9	PS-10
	9	AR-516-9							PS-7	PS-8	PS-9	PS-10
	10	AR-516-10								PS-8	PS-9	PS-10
CR Poly w/CR Femur	1	AR-517-1	CR-1	CR-2	CR-3							
	2	AR-517-2	CR-1	CR-2	CR-3	CR-4						
	3	AR-517-3	CR-1	CR-2	CR-3	CR-4	CR-5					
	4	AR-517-4		CR-2	CR-3	CR-4	CR-5	CR-6				
	5	AR-517-5			CR-3	CR-4	CR-5	CR-6	CR-7			
	6	AR-517-6				CR-4	CR-5	CR-6	CR-7	CR-8		
	7	AR-517-7					CR-5	CR-6	CR-7	CR-8	CR-9	
	8	AR-517-8						CR-6	CR-7	CR-8	CR-9	CR-10
	9	AR-517-9							CR-7	CR-8	CR-9	CR-10
	10	AR-517-10								CR-8	CR-9	CR-10
CR Plus Poly w/CR Femur	1	AR-517-1	CR Plus-1	CR Plus-2								
	2	AR-517-2	CR Plus-1	CR Plus-2	CR Plus-3							
	3	AR-517-3		CR Plus-2	CR Plus-3	CR Plus-4						
	4	AR-517-4			CR Plus-3	CR Plus-4	CR Plus-5					
	5	AR-517-5				CR Plus-4	CR Plus-5	CR Plus-6				
	6	AR-517-6					CR Plus-5	CR Plus-6	CR Plus-7			
	7	AR-517-7						CR Plus-6	CR Plus-7	CR Plus-8		
	8	AR-517-8							CR Plus-7	CR Plus-8	CR Plus-9	
	9	AR-517-9								CR Plus-8	CR Plus-9	CR Plus-10
	10	AR-517-10									CR Plus-9	CR Plus-10

- PS Bearing (sizes 8 mm–20 mm)
- CR Bearing (sizes 8 mm–18 mm)
- CR Plus Bearing (sizes 8 mm–18 mm)



# iBalance<sup>®</sup> TKA

## Implant Sizing Compatibility Matrix

