


JUMPSTART® ANTIMICROBIAL WOUND DRESSING


JumpStart

ADVANCED MICROCURRENT TECHNOLOGY®

New Knee Dressing




Shorter overall dressing length for improved fit and mobility.



New Energel™ Wound Hydrogel

- Compatible with JumpStart antimicrobial wound dressing
- Maintains JumpStart wound dressing's electric field for 5- to 7-day wear time
- Creates a moist wound environment
- Double packaged sterile sachet for use in operating room
- Nonwasteful 7.5 g sachet for true single use
- Available in a box of 10 sachets



Arthrex®
www.Arthrex.com

To order, call 1.800.934.4404

Advanced Microcurrent Technology® and Energel™ are registered trademarks of Vomar Innovations, Inc.

*Utility and design patents pending

© 2018 Arthrex, Inc. All rights reserved. LS1-00069-EN_A