


# TPLO Plate System

## Weight Reference Chart



InternalBrace™ Anchors	TPLO Plate Sizes						InternalBrace Sutures
2.4 mm/2.5 mm PushLock® anchor	2.0 mm						SutureTape suture, 1.3 mm
2.5 mm/2.9 mm PushLock anchor		2.4 mm					SutureTape suture, 1.3 mm
3.5 mm SwiveLock® anchor			2.7 mm				LabralTape™ suture, 1.5 mm FiberTape® suture, 2.0 mm
3.5 mm/4.75 mm SwiveLock anchor				3.5 mm, small			LabralTape suture, 1.5 mm FiberTape suture, 2.0 mm
4.75 mm SwiveLock anchor					3.5 mm, standard		FiberTape suture, 2.0 mm
4.75 mm/ 5.5 mm SwiveLock anchor						3.5 mm, broad	FiberTape suture, 2.0 mm
5.5 mm SwiveLock anchor							4.5 mm FiberTape suture, 2.0 mm

All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultant's recommendation. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.